



# Orthopedic Trauma Solutions



[www.4umedical.com.tr](http://www.4umedical.com.tr)



**Medical**

Keeps You  
Moving



**PLEASE SCAN QR CODE  
TO LEARN MORE!**



# Contents

## 1. Trauma Locking Plate System

- Upper Limb Plates.....
- Clavicula Shaft Locking Plate.....
- Clavicula Distal Locking Plate.....
- Clavicula Hook Locking Plate.....
- Humerus Proximal Locking Plate.....
- Humerus Distal Lateral Locking Plate.....
- Humerus Distal Posterolateral Locking Plate.....
- Humerus Distal Medial Locking Plate.....
- Olecranon Locking Plate.....
- Olecranon Hook Locking Plate.....
- Small 3.5 Compression Plate.....
- Small 3.5 Tp2 Locking Straight Plate.....
- Tubular Locking Plate.....
- Proximal Radius Locking Plate.....
- Distal Radius Volar Locking Plate.....
- Distal Radius Volar Locking Plate Small.....
- Distal Radius Volar Locking Plate Tp2.....
- Radius Distal Dorsal T Locking Plate.....
- Radius Distal Dorsal L Locking Plate.....
- Distal Radius Styloid Plate.....
- Radius Distal/Ulna Hook Locking Plate.....
- Ulna Proximal Hook Locking Plate.....
- Distal Ulna Locking Plate.....
- Coronoid Locking Plate.....
- Metacarpal Hand/Wrist Arthrodesis Plate.....

## 1. Trauma Locking Plate System

- Lower Limb Plates.....
- Tubular Locking Plate.....
- Femur Proximal Locking Plate.....
- Femur Distal Locking Plate.....
- Femur High Osteotomy Plate 4.5/5.0 mm.....
- Tibia Distal Medial Locking Plate.....
- Tibia Distal Medial Locking Plate 3.5.....
- Large Broad Locking Plate.....
- Large Broad Locking T Plate.....
- Large Broad Locking T Plate.....
- Large Broad T-Support Locking Plate.....
- Tibia High Osteotomy Plate 4.5/5.0 mm.....
- Large Broad L-Support Locking Plate.....
- Proximal Medial/Anteromedial Tibia Locking Plate.....
- Tibia Proximal Lateral Locking Plate.....
- Tibia Distal Lateral Locking Plate.....
- Tibia Distal Lateral Locking Plate Tp2.....
- Pilon Locking Plate.....
- Fibula Locking Plate.....
- Reconstruction Shaft Locking Plate.....
- Calcaneus Locking Plate.....
- Fibula Distal Locking Plate.....
- Small Femur Proximal Hip Locking Plate.....
- Pelvic Reconstruction Curved Locking Plate.....
- Pelvic Symphyseal Plate.....
- Pelvic Reconstruction Straight Locking Plate.....
- Eight Non Locking Plate.....





# Contents

## 2. Trauma Screw

- Locking Self Tapping Cortical Screw.....
  - Locking Self Tapping Cortical Screw Ø 1,5 mm.....
  - Locking Self Tapping Cortical Screw Ø 2,0 mm.....
  - Locking Self Tapping Cortical Screw Ø 2,5 mm.....
  - Locking Self Tapping Cortical Screw Ø 2,7 mm.....
  - Locking Self Tapping Cortical Screw Ø 3,5 mm.....
  - Locking Self Tapping Cortical Screw Ø 5,0 mm.....
- Cortical Self Tapping Screw.....
  - Cortical Self Tapping Screw Ø 1,5 mm.....
  - Cortical Self Tapping Screw Ø 2,0 mm.....
  - Cortical Self Tapping Screw Ø 2,5 mm.....
  - Cortical Self Tapping Screw Ø 2,7 mm.....
  - Cortical Self Tapping Screw Ø 3,5 mm.....
  - Cortical Self Tapping Screw Ø 4,5 mm.....
- Locking Cancellous Spongiose Screw.....
  - Locking Cancellous/Spongiose Screw Ø 4,0 mm.....
  - Locking Cancellous/Spongiose Screw Ø 5,0 mm.....
  - Cancellous/Spongiose Screw Full Thread Ø 4,0 mm.....
  - Cancellous/Spongiose Screw Full Thread Ø 6,5 mm.....
  - Cancellous Screw Half Thread Ø 4,0 mm.....
  - Cancellous Screw 16 mm Half Thread Ø 6,5 mm.....
  - Cancellous Screw 32 mm Thread Ø 6,5 mm.....

## 2. Trauma Screw

- Cannulated Compression Screw.....
  - Cannulated Compression Screw Ø 2,7 mm.....
  - Cannulated Compression Screw Full Thread Ø 2,7 mm.....
  - Cannulated Compression Screw Ø 3,5 mm.....
  - Cannulated Compression Screw Ø 4,5 mm.....
  - Cannulated Compression Screw 16 mm Thread Ø 6,5 mm.....
  - Cannulated Compression Screw 32 mm Thread Ø 6,5 mm.....
  - Cannulated Compression Screw 16 mm Thread Ø 7,3 mm.....
  - Cannulated Compression Screw 32 mm Thread Ø 7,3 mm.....
  - Locking Cannulated Full Thread Screw Ø 2,5 mm.....
  - Locking Cannulated Full Thread Screw Ø 3,5 mm.....
  - Locking Cannulated Full Thread Screw Ø 5,0 mm.....
  - Locking Cannulated Full Thread Screw Ø 6,5 mm.....
  - Washer.....
- Headless Cannulated Screw.....
  - Headless Cannulated Full Thread Screw Ø 2,5 mm.....
  - Headless Cannulated Full Thread Screw Ø 3,5 mm.....
  - Headless Cannulated Full Thread Screw Ø 4,0 mm.....
  - Headless Cannulated Full Thread Screw Ø 5,0 mm.....
- Herbert Screw.....
  - Herbert Screw Ø 2,5 mm.....
  - Herbert Screw Ø 3,0 mm.....
  - Herbert Screw Ø 5,0 mm.....



# Contents

|                                                           |  |
|-----------------------------------------------------------|--|
| • Cable Locking Screw.....                                |  |
| • Cap Screw 2,5/3,5/5,0.....                              |  |
| 3. Mini Plate System.....                                 |  |
| • 1,5 mm Hand / Foot Mini Locking Plate Series.....       |  |
| • 1,5 mm Locking / Non - Locking Screw.....               |  |
| • 2,0 mm Hand / Foot Mini Locking Plate Series.....       |  |
| • 2,0 mm Locking / Non - Locking Screw.....               |  |
| • 2,5 mm Hand / Foot Mini Locking Plate Series.....       |  |
| • 2,5 mm Locking / Non - Locking Screw.....               |  |
| • 2,7 mm Hand / Foot Mini Locking Plate Series.....       |  |
| • 2,7 mm Hand / Foot Mini Locking Plate Series.....       |  |
| • Metacarpal / Metatarsal Fusion Plate Series.....        |  |
| • 2,5 mm Locking / Non - Locking Screw.....               |  |
| 4. Nail System.....                                       |  |
| • Humerus Cannulated Nail / Nail Screw.....               |  |
| • PFNA Nail.....                                          |  |
| • PFNA Nail Screw.....                                    |  |
| • Intramedullary Cannulated Lateral Femur Nail.....       |  |
| • Intramedullary Cannulated Lateral Femur Nail Screw..... |  |
| • Tibial Nail Cannulated.....                             |  |
| • Tibial Nail Cannulated Screw.....                       |  |
| • Arthrodesis Nail / Screw.....                           |  |
| • Flexible Intramedullary Nail / End Cap.....             |  |



# Introduction



4U Medical Ltd., Co. aims a pioneering leader in sales of high-quality orthopedic trauma solutions. Specializing in advanced implants and instruments, our mission is to enhance patient outcomes and support medical professionals in the field of orthopedic surgery.

Our comprehensive product portfolio includes state-of-the-art devices for osteosynthesis, with a focus on the treatment of fractures and other trauma-related conditions. From intramedullary nails to dynamic compression systems, our products are engineered to ensure precision, durability, and optimal healing.


Key offerings include:

- Trauma Screws: A full range of screws (cortical, cancellous, locking, and non-locking) designed to meet the specific needs of fracture stabilization.
- Limited Contact Plates: Innovative plates designed to minimize soft tissue irritation while ensuring secure fixation.
- Fracture Reduction Tools: Specialized instruments to aid in the realignment and stabilization of fractured bones.
- Osteosynthesis Systems: Comprehensive solutions for bone repair and regeneration, supporting complex surgical procedures.

By combining cutting-edge technology with rigorous quality control, 4U Medical Ltd., Co. aims to provide surgeons with reliable tools to restore mobility and improve the quality of life for patients. Partner with us for innovative solutions in orthopedic trauma care.

4U Medical Ltd., Co. – Keeps You Moving.

## CLAVICULA SHAFT LOCKING PLATE



**R**

Thickness: 3  
Width: 10

### Screws

Ø3.5 mm – Locking Self Tapping Screw  
Ø3.5 mm – Cortical Self Tapping Screw  
Ø4,0 mm – Locking Cancellous Screw  
Ø4.0 mm – Cancellous Screw

### Code Holes Side Length

|              |    |   |     |
|--------------|----|---|-----|
| 4U1011213504 | 4  | R | 49  |
| 4U1011213505 | 5  | R | 60  |
| 4U1011213506 | 6  | R | 72  |
| 4U1011213507 | 7  | R | 84  |
| 4U1011213508 | 8  | R | 96  |
| 4U1011213509 | 9  | R | 107 |
| 4U1011213510 | 10 | R | 118 |
| 4U1011213511 | 11 | R | 129 |
| 4U1011213512 | 12 | R | 140 |
| 4U1011223504 | 4  | L | 49  |
| 4U1011223505 | 5  | L | 60  |
| 4U1011223506 | 6  | L | 72  |
| 4U1011223507 | 7  | L | 84  |
| 4U1011223508 | 8  | L | 96  |
| 4U1011223509 | 9  | L | 107 |
| 4U1011223510 | 10 | L | 118 |
| 4U1011223511 | 11 | L | 129 |
| 4U1011223512 | 12 | L | 140 |

## CLAVICULA DISTAL LOCKING PLATE

### Code Holes Side Length

|              |    |   |     |
|--------------|----|---|-----|
| 4U1011313503 | 3  | R | 54  |
| 4U1011313504 | 4  | R | 67  |
| 4U1011313505 | 5  | R | 80  |
| 4U1011313506 | 6  | R | 93  |
| 4U1011313507 | 7  | R | 106 |
| 4U1011313508 | 8  | R | 119 |
| 4U1011313509 | 9  | R | 132 |
| 4U1011313510 | 10 | R | 145 |
| 4U1011323503 | 3  | L | 54  |
| 4U1011323504 | 4  | L | 67  |
| 4U1011323505 | 5  | L | 80  |
| 4U1011323506 | 6  | L | 93  |
| 4U1011323507 | 7  | L | 106 |
| 4U1011323508 | 8  | L | 119 |
| 4U1011323509 | 9  | L | 132 |
| 4U1011323510 | 10 | L | 145 |

### Screws

Ø 2.5 mm – Locking Self Tapping Screw  
Ø 3.5 mm – Locking Self Tapping Screw  
Ø 3.5 mm – Cortical Self Tapping Screw  
Ø 4,0 mm – Locking Cancellous Screw  
Ø 4.0 mm – Cancellous Screw



**R**

Thickness: 3,4  
Width: 10

## CLAVICULA HOOK LOCKING PLATE



**R**

Thickness: 3,5  
Width: 11

### Screws

Ø 3.5 mm – Locking Self Tapping Screw  
Ø 3.5 mm – Cortical Self Tapping Screw  
Ø 4,0 mm – Locking Cancellous Screw  
Ø 4.0 mm – Cancellous Screw

### Code Holes Side Length

|              |   |   |     |
|--------------|---|---|-----|
| 4U1011513503 | 3 | R | 57  |
| 4U1011513504 | 4 | R | 72  |
| 4U1011513505 | 5 | R | 87  |
| 4U1011513506 | 6 | R | 102 |
| 4U1011513507 | 7 | R | 117 |
| 4U1011513508 | 8 | R | 132 |
| 4U1011523503 | 3 | L | 57  |
| 4U1011523504 | 4 | L | 72  |
| 4U1011523505 | 5 | L | 87  |
| 4U1011523506 | 6 | L | 102 |
| 4U1011523507 | 7 | L | 117 |
| 4U1011523508 | 8 | L | 132 |

# HUMERUS PROXIMAL LOCKING PLATE

Thickness: 4  
Width: 12



### Screws

- Ø 3.5 mm – Locking Self Tapping Screw
- Ø 3.5 mm – Cortical Self Tapping Screw
- Ø 4,0 mm – Locking Cancellous Screw
- Ø 4.0 mm – Cancellous Screw

| Code | Holes | Side | Length |
|------|-------|------|--------|
|------|-------|------|--------|

|              |    |   |     |
|--------------|----|---|-----|
| 4U1021103503 | 3  | - | 105 |
| 4U1021103504 | 4  | - | 122 |
| 4U1021103505 | 5  | - | 139 |
| 4U1021103506 | 6  | - | 156 |
| 4U1021103507 | 7  | - | 173 |
| 4U1021103508 | 8  | - | 190 |
| 4U1021103509 | 9  | - | 207 |
| 4U1021103510 | 10 | - | 224 |
| 4U1021103511 | 11 | - | 241 |
| 4U1021103512 | 12 | - | 258 |
| 4U1021103513 | 13 | - | 275 |

# HUMERUS DISTAL LATERAL LOCKING PLATE

| Code | Holes | Side | Length |
|------|-------|------|--------|
|------|-------|------|--------|

|              |    |   |     |
|--------------|----|---|-----|
| 4U1021313503 | 3  | R | 85  |
| 4U1021313504 | 4  | R | 98  |
| 4U1021313505 | 5  | R | 111 |
| 4U1021313506 | 6  | R | 124 |
| 4U1021313507 | 7  | R | 137 |
| 4U1021313508 | 8  | R | 150 |
| 4U1021313509 | 9  | R | 163 |
| 4U1021313510 | 10 | R | 176 |
| 4U1021313511 | 11 | R | 189 |
| 4U1021313512 | 12 | R | 202 |
| 4U1021313513 | 13 | R | 215 |
| 4U1021323503 | 3  | L | 85  |
| 4U1021323504 | 4  | L | 98  |
| 4U1021323505 | 5  | L | 111 |
| 4U1021323506 | 6  | L | 124 |
| 4U1021323507 | 7  | L | 137 |
| 4U1021323508 | 8  | L | 150 |
| 4U1021323509 | 9  | L | 163 |
| 4U1021323510 | 10 | L | 176 |
| 4U1021323511 | 11 | L | 189 |
| 4U1021323512 | 12 | L | 202 |
| 4U1021323513 | 13 | L | 215 |

### Screws

- Ø3.5 mm – Locking Self Tapping Screw
- Ø3.5 mm – Cortical Self Tapping Screw
- Ø4,0 mm – Locking Cancellous Screw
- Ø4.0 mm – Cancellous Screw

R  
Thickness: 3,5  
Width: 11



# HUMERUS DISTAL POSTEROLATERAL LOCKING PLATE

R  
Thickness: 3  
Width: 12




### Screws

- Ø 2.5 mm – Locking Self Tapping Screw
- Ø 2.5 mm – Cortical Self Tapping Screw
- Ø 3.5 mm – Locking Self Tapping Screw
- Ø 3.5 mm – Cortical Self Tapping Screw
- Ø 4,0 mm – Locking Cancellous Screw
- Ø 4.0 mm – Cancellous Screw

| Code | Holes | Side | Length |
|------|-------|------|--------|
|------|-------|------|--------|

|              |    |   |     |
|--------------|----|---|-----|
| 4U1021613503 | 3  | R | 74  |
| 4U1021613504 | 4  | R | 85  |
| 4U1021613505 | 5  | R | 98  |
| 4U1021613506 | 6  | R | 111 |
| 4U1021613507 | 7  | R | 124 |
| 4U1021613508 | 8  | R | 137 |
| 4U1021613509 | 9  | R | 150 |
| 4U1021613510 | 10 | R | 163 |
| 4U1021613511 | 11 | R | 176 |
| 4U1021613512 | 12 | R | 189 |
| 4U1021613513 | 13 | R | 202 |
| 4U1021623503 | 3  | L | 74  |
| 4U1021623504 | 4  | L | 85  |
| 4U1021623505 | 5  | L | 98  |
| 4U1021623506 | 6  | L | 111 |
| 4U1021623507 | 7  | L | 124 |
| 4U1021623508 | 8  | L | 137 |
| 4U1021623509 | 9  | L | 150 |
| 4U1021623510 | 10 | L | 163 |
| 4U1021623511 | 11 | L | 176 |
| 4U1021623512 | 12 | L | 189 |
| 4U1021623513 | 13 | L | 202 |

## HUMERUS DISTAL MEDIAL LOCKING PLATE



**R**

Thickness: 3  
Width: 11

### Screws

- Ø 2,5 mm – Locking Self Tapping Screw
- Ø 2,5 mm – Cortical Self Tapping Screw
- Ø 3.5 mm – Locking Self Tapping Screw
- Ø 3.5 mm – Cortical Self Tapping Screw
- Ø 4,0 mm – Locking Cancellous Screw
- Ø 4.0 mm – Cancellous Screw

| Code         | Holes | Side | Length |
|--------------|-------|------|--------|
| 4U1021713503 | 3     | R    | 70     |
| 4U1021713504 | 4     | R    | 82     |
| 4U1021713505 | 5     | R    | 95     |
| 4U1021713506 | 6     | R    | 108    |
| 4U1021713507 | 7     | R    | 121    |
| 4U1021713508 | 8     | R    | 134    |
| 4U1021713509 | 9     | R    | 147    |
| 4U1021713510 | 10    | R    | 160    |
| 4U1021713511 | 11    | R    | 173    |
| 4U1021713512 | 12    | R    | 186    |
| 4U1021713513 | 13    | R    | 199    |
| 4U1021723503 | 3     | L    | 70     |
| 4U1021723504 | 4     | L    | 82     |
| 4U1021723505 | 5     | L    | 95     |
| 4U1021723506 | 6     | L    | 108    |
| 4U1021723507 | 7     | L    | 121    |
| 4U1021723508 | 8     | L    | 134    |
| 4U1021723509 | 9     | L    | 147    |
| 4U1021723510 | 10    | L    | 160    |
| 4U1021723511 | 11    | L    | 173    |
| 4U1021723512 | 12    | L    | 186    |
| 4U1021723513 | 13    | L    | 199    |

## OLECRANON LOCKING PLATE

### Code Holes Side Length

|              |    |   |     |
|--------------|----|---|-----|
| 4U1031113503 | 3  | R | 101 |
| 4U1031113504 | 4  | R | 110 |
| 4U1031113505 | 5  | R | 125 |
| 4U1031113506 | 6  | R | 138 |
| 4U1031113507 | 7  | R | 151 |
| 4U1031113508 | 8  | R | 164 |
| 4U1031113509 | 9  | R | 177 |
| 4U1031113510 | 10 | R | 190 |
| 4U1031113511 | 11 | R | 203 |
| 4U1031113512 | 12 | R | 216 |
| 4U1031113513 | 13 | R | 229 |
| 4U1031123503 | 3  | L | 101 |
| 4U1031123504 | 4  | L | 110 |
| 4U1031123505 | 5  | L | 125 |
| 4U1031123506 | 6  | L | 138 |
| 4U1031123507 | 7  | L | 151 |
| 4U1031123508 | 8  | L | 164 |
| 4U1031123509 | 9  | L | 177 |
| 4U1031123510 | 10 | L | 190 |
| 4U1031123511 | 11 | L | 203 |
| 4U1031123512 | 12 | L | 216 |
| 4U1031123513 | 13 | L | 229 |

### Screws

- Ø 3.5 mm – Locking Self Tapping Screw
- Ø 3.5 mm – Cortical Self Tapping Screw
- Ø 4,0 mm – Locking Cancellous Screw
- Ø 4.0 mm – Cancellous Screw



**R**

Thickness: 3  
Width: 10

## OLECRANON HOOK LOCKING PLATE



Thickness: 3  
Width: 10

### Screws


- Ø 3.5 mm – Locking Self Tapping Screw
- Ø 3.5 mm – Cortical Self Tapping Screw
- Ø 4,0 mm – Locking Cancellous Screw
- Ø 4.0 mm – Cancellous Screw

### Code Holes Side Length

|              |   |   |     |
|--------------|---|---|-----|
| 4U8081203503 | 3 | - | 53  |
| 4U8081203504 | 4 | - | 67  |
| 4U8081203505 | 5 | - | 81  |
| 4U8081203506 | 6 | - | 95  |
| 4U8081203507 | 7 | - | 109 |
| 4U8081203508 | 8 | - | 123 |
| 4U8081203509 | 9 | - | 137 |

## SMALL 3.5 COMPRESSION PLATE

Thickness: 3,5  
Width: 12



### Screws

Ø 3.5 mm – Locking Self Tapping Screw  
Ø 3.5 mm – Cortical Self Tapping Screw  
Ø 4,0 mm – Locking Cancellous Screw  
Ø 4.0 mm – Cancellous Screw

| Code         | Holes | Side | Length |
|--------------|-------|------|--------|
| 4U1031283524 | 4     | -    | 62     |
| 4U1031283515 | 5     | -    | 74,5   |
| 4U1031283516 | 6     | -    | 87     |
| 4U1031283517 | 7     | -    | 99,5   |
| 4U1031283518 | 8     | -    | 112    |
| 4U1031283519 | 9     | -    | 124,5  |
| 4U1031283510 | 10    | -    | 137    |
| 4U1031283511 | 11    | -    | 149,5  |
| 4U1031283512 | 12    | -    | 162    |
| 4U1031283513 | 13    | -    | 174,5  |
| 4U1031283514 | 14    | -    | 187    |

## SMALL 3.5 TP2 LOCKING STRAIGHT PLATE

### Code Holes Side Length

|              |    |   |       |
|--------------|----|---|-------|
| 4U1031383504 | 4  | - | 62    |
| 4U1031383505 | 5  | - | 74,5  |
| 4U1031383506 | 6  | - | 87    |
| 4U1031383507 | 7  | - | 99,5  |
| 4U1031383508 | 8  | - | 112   |
| 4U1031383509 | 9  | - | 124,5 |
| 4U1031383510 | 10 | - | 137   |
| 4U1031383511 | 11 | - | 149,5 |
| 4U1031383512 | 12 | - | 162   |
| 4U1031383513 | 13 | - | 174,5 |
| 4U1031383514 | 14 | - | 187   |

### Screws

Ø 3.5 mm – Locking Self Tapping Screw  
Ø 3.5 mm – Cortical Self Tapping Screw  
Ø 4,0 mm – Locking Cancellous Screw  
4.0 mm – Cancellous Screw



Thickness: 3,5  
Width: 13

## TUBULAR LOCKING PLATE

Thickness: 1,4  
Width: 10



### Screws

Ø 3.5 mm – Locking Self Tapping Screw  
Ø 3.5 mm – Cortical Self Tapping Screw  
Ø 4,0 mm – Locking Cancellous Screw  
4.0 mm – Cancellous Screw

### Code Holes Side Length

|              |    |   |     |
|--------------|----|---|-----|
| 4U1081203504 | 4  | - | 54  |
| 4U1081203505 | 5  | - | 66  |
| 4U1081203506 | 6  | - | 78  |
| 4U1081203507 | 7  | - | 90  |
| 4U1081203508 | 8  | - | 102 |
| 4U1081203509 | 9  | - | 114 |
| 4U1081203510 | 10 | - | 126 |
| 4U1081203511 | 11 | - | 138 |
| 4U1081203512 | 12 | - | 150 |

# PROXIMAL RADIUS LOCKING PLATE



**R**  
Thickness: 1,9  
Width: 9

## Screws

Ø 2.5 mm – Locking Self Tapping Screw

Ø 2.5 mm – Cortical Self Tapping Screw

| Code         | Holes | Side | Length |
|--------------|-------|------|--------|
| 4U1031112502 | 2     | R    | 33     |
| 4U1031112503 | 3     | R    | 42     |
| 4U1031112504 | 4     | R    | 51     |
| 4U1031112505 | 5     | R    | 60     |
| 4U1031112506 | 6     | R    | 69     |
| 4U1031122502 | 2     | L    | 33     |
| 4U1031122503 | 3     | L    | 42     |
| 4U1031122504 | 4     | L    | 51     |
| 4U1031122505 | 5     | L    | 60     |
| 4U1031122506 | 6     | L    | 69     |

# DISTAL RADIUS VOLAR LOCKING PLATE

| Code         | Holes | Side | Length |
|--------------|-------|------|--------|
| 4U1031312502 | 2     | R    | 51     |
| 4U1031312503 | 3     | R    | 59     |
| 4U1031312504 | 4     | R    | 67     |
| 4U1031312505 | 5     | R    | 75     |
| 4U1031312506 | 6     | R    | 83     |
| 4U1031312507 | 7     | R    | 91     |
| 4U1031312508 | 8     | R    | 99     |
| 4U1031312509 | 9     | R    | 107    |
| 4U1031322502 | 2     | L    | 51     |
| 4U1031322503 | 3     | L    | 59     |
| 4U1031322504 | 4     | L    | 67     |
| 4U1031322505 | 5     | L    | 75     |
| 4U1031322506 | 6     | L    | 83     |
| 4U1031322507 | 7     | L    | 91     |
| 4U1031322508 | 8     | L    | 99     |
| 4U1031322509 | 9     | L    | 107    |

## Screws

Ø 2.5 mm – Locking Self Tapping Screw

Ø 2.5 mm – Cortical Self Tapping Screw



**R**  
Thickness: 2  
Width: 8,5

# DISTAL RADIUS VOLAR LOCKING PLATE SMALL



**R**  
Thickness: 2  
Width: 8,5

## Screws

Ø 2.5 mm – Locking Self Tapping Screw

Ø 2.5 mm – Cortical Self Tapping Screw

| Code         | Holes | Side | Length |
|--------------|-------|------|--------|
| 4U1031332502 | 2     | R    | 51     |
| 4U1031332503 | 3     | R    | 59     |
| 4U1031332504 | 4     | R    | 67     |
| 4U1031332505 | 5     | R    | 75     |
| 4U1031332506 | 6     | R    | 83     |
| 4U1031332507 | 7     | R    | 91     |
| 4U1031332508 | 8     | R    | 99     |
| 4U1031332509 | 9     | R    | 107    |
| 4U1031342502 | 2     | L    | 51     |
| 4U1031342503 | 3     | L    | 59     |
| 4U1031342504 | 4     | L    | 67     |
| 4U1031342505 | 5     | L    | 75     |
| 4U1031342506 | 6     | L    | 83     |
| 4U1031342507 | 7     | L    | 91     |
| 4U1031342508 | 8     | L    | 99     |
| 4U1031342509 | 9     | L    | 107    |

## DISTAL RADIUS VOLAR LOCKING PLATE TP2



**R**

Thickness: 2  
Width: 8,5

Screws

Ø 2.5 mm – Locking Self Tapping Screw

Ø 2.5 mm – Cortical Self Tapping Screw

| Code | Holes | Side | Length |
|------|-------|------|--------|
|------|-------|------|--------|

|              |   |   |     |
|--------------|---|---|-----|
| 4U1031352502 | 2 | R | 51  |
| 4U1031352503 | 3 | R | 59  |
| 4U1031352504 | 4 | R | 67  |
| 4U1031352505 | 5 | R | 75  |
| 4U1031352506 | 6 | R | 83  |
| 4U1031352507 | 7 | R | 91  |
| 4U1031352508 | 8 | R | 99  |
| 4U1031352509 | 9 | R | 107 |
| 4U1031362502 | 2 | L | 51  |
| 4U1031362503 | 3 | L | 59  |
| 4U1031362504 | 4 | L | 67  |
| 4U1031362505 | 5 | L | 75  |
| 4U1031362506 | 6 | L | 83  |
| 4U1031362507 | 7 | L | 91  |
| 4U1031362508 | 8 | L | 99  |
| 4U1031362509 | 9 | L | 107 |

## RADIUS DISTAL DORSAL T LOCKING PLATE

| Code | Holes | Side | Length |
|------|-------|------|--------|
|------|-------|------|--------|

|              |   |   |    |
|--------------|---|---|----|
| 4U1031402503 | 3 | - | 34 |
| 4U1031402504 | 4 | - | 40 |
| 4U1031402505 | 5 | - | 46 |
| 4U1031402506 | 6 | - | 52 |



Thickness: 2  
Width: 7

Screws

Ø 2.5 mm – Locking Self Tapping Screw

Ø 2.5 mm – Cortical Self Tapping Screw

## RADIUS DISTAL DORSAL L LOCKING PLATE

**R**

Thickness: 1,8  
Width: 7,5



Screws

Ø 2.5 mm – Locking Self Tapping Screw

Ø 2.5 mm – Cortical Self Tapping Screw

| Code | Holes | Side | Length |
|------|-------|------|--------|
|------|-------|------|--------|

|              |   |   |    |
|--------------|---|---|----|
| 4U1031412505 | 3 | R | 42 |
| 4U1031412506 | 4 | R | 48 |
| 4U1031412507 | 5 | R | 54 |
| 4U1031412508 | 6 | R | 60 |
| 4U1031422501 | 3 | L | 42 |
| 4U1031422502 | 4 | L | 48 |
| 4U1031422503 | 5 | L | 54 |
| 4U1031422504 | 6 | L | 60 |

## DISTAL RADIUS STYLOID PLATE

| Code | Holes | Side | Length |
|------|-------|------|--------|
|------|-------|------|--------|

|              |   |   |      |
|--------------|---|---|------|
| 4U1031402511 | 4 | - | 68,5 |
|--------------|---|---|------|



Thickness: 1,8  
Width: 7,5

Screws

Ø 2.5 mm – Locking Self Tapping Screw

Ø 2.5 mm – Cortical Self Tapping Screw

# RADIUS DISTAL/ULNA HOOK LOCKING PLATE

Thickness: 1,7  
Width: 6,2



### Screws

- Ø 2.5 mm – Locking Self Tapping Screw
- Ø 2.5 mm – Cortical Self Tapping Screw

| Code         | Holes | Side | Length |
|--------------|-------|------|--------|
| 4U1031402510 | 6     | -    | 42     |

# ULNA PROXIMAL HOOK LOCKING PLATE

## Code Holes Side Length

|              |    |   |     |
|--------------|----|---|-----|
| 4U1041103503 | 3  | - | 53  |
| 4U1041103504 | 4  | - | 67  |
| 4U1041103505 | 5  | - | 81  |
| 4U1041103506 | 6  | - | 95  |
| 4U1041103507 | 7  | - | 109 |
| 4U1041103508 | 8  | - | 123 |
| 4U1041103509 | 9  | - | 137 |
| 4U1041103510 | 10 | - | 151 |

### Screws

- Ø 3.5 mm – Locking Self Tapping Screw
- Ø 3.5 mm – Cortical Self Tapping Screw
- Ø 4,0 mm – Locking Cancellous Screw
- Ø 4.0 mm – Cancellous Screw

Thickness: 3  
Width: 10



# DISTAL ULNA LOCKING PLATE

**R**  
Thickness: 2,5  
Width: 9



### Screws

- Ø 2.5 mm – Locking Self Tapping Screw
- Ø 2.7 mm – Locking Self Tapping Screw
- Ø 3.5 mm – Cortical Self Tapping Screw
- Ø 3.5 mm – Locking Self Tapping Screw

## Code Holes Side Length

|              |   |   |     |
|--------------|---|---|-----|
| 4U1041313503 | 3 | R | 49  |
| 4U1041313504 | 4 | R | 60  |
| 4U1041313505 | 5 | R | 71  |
| 4U1041313506 | 6 | R | 86  |
| 4U1041313507 | 7 | R | 97  |
| 4U1041313508 | 8 | R | 112 |
| 4U1041313509 | 9 | R | 123 |
| 4U1041323503 | 3 | L | 49  |
| 4U1041323504 | 4 | L | 60  |
| 4U1041323505 | 5 | L | 71  |
| 4U1041323506 | 6 | L | 86  |
| 4U1041323507 | 7 | L | 97  |
| 4U1041323508 | 8 | L | 112 |
| 4U1041323509 | 9 | L | 123 |

# CORONOID LOCKING PLATE

## Code Holes Side Length

|              |   |   |        |
|--------------|---|---|--------|
| 4U1001000150 | 2 | R | Small  |
| 4U1001000151 | 2 | R | Medium |
| 4U1001000152 | 2 | R | Large  |
| 4U1001000250 | 2 | L | Small  |
| 4U1001000251 | 2 | L | Medium |
| 4U1001000252 | 2 | L | Large  |

### Screws

- Ø 3.5 mm – Locking Self Tapping Screw
- Ø 3.5 mm – Cortical Self Tapping Screw
- Ø 4,0 mm – Locking Cancellous Screw
- Ø 4.0 mm – Cancellous Screw

**R**  
Thickness: 2,0  
Width: 9,0



# METACARPAL HAND/WRIST ARTHRODESIS PLATE

Code

4U1031902500

Size

S



Thickness: 3  
Width: 10

## Screws

- Ø 2.5 mm – Locking Self Tapping Screw
- Ø 2.7 mm – Locking Self Tapping Screw
- Ø 3.5 mm – Locking Self Tapping Screw
- Ø 3.5 mm – Cortical Self Tapping Screw
- Ø 4.0 mm – Cancellous Screw

Code


4U1031902501  
4U1031902502

Size

M  
L



# FEMUR PROXIMAL LOCKING PLATE



| Code         | Holes | Side | Length |
|--------------|-------|------|--------|
| 4U2061115503 | 3     | R    | 120    |
| 4U2061115504 | 4     | R    | 138    |
| 4U2061115505 | 5     | R    | 156    |
| 4U2061115506 | 6     | R    | 174    |
| 4U2061115507 | 7     | R    | 192    |
| 4U2061115508 | 8     | R    | 210    |
| 4U2061115509 | 9     | R    | 228    |
| 4U2061115510 | 10    | R    | 246    |
| 4U2061115511 | 11    | R    | 264    |
| 4U2061115512 | 12    | R    | 282    |
| 4U2061115513 | 13    | R    | 300    |
| 4U2061127003 | 3     | L    | 120    |
| 4U2061127004 | 4     | L    | 138    |
| 4U2061127005 | 5     | L    | 156    |
| 4U2061127006 | 6     | L    | 174    |
| 4U2061127007 | 7     | L    | 192    |
| 4U2061127008 | 8     | L    | 210    |
| 4U2061127009 | 9     | L    | 228    |
| 4U2061127010 | 10    | L    | 246    |
| 4U2061127011 | 11    | L    | 264    |
| 4U2061127012 | 12    | L    | 282    |
| 4U2061127013 | 13    | L    | 300    |


Screws  
Ø 5.0 mm – Locking Self Tapping Screw Ø  
4.5 mm – Cortical Self Tapping Screw Ø 5.0  
mm – Locking Cancellous Screw Ø 6.5 mm –  
Locking Cannulated Full Thread Screw

# FEMUR DISTAL LOCKING PLATE

| Code         | Holes | Side | Length |
|--------------|-------|------|--------|
| 4U2061316005 | 5     | R    | 160    |
| 4U2061316006 | 6     | R    | 178    |
| 4U2061316007 | 7     | R    | 196    |
| 4U2061316008 | 8     | R    | 214    |
| 4U2061316009 | 9     | R    | 232    |
| 4U2061316010 | 10    | R    | 250    |
| 4U2061316011 | 11    | R    | 268    |
| 4U2061316012 | 12    | R    | 286    |
| 4U2061316013 | 13    | R    | 304    |
| 4U2061325005 | 5     | L    | 160    |
| 4U2061325006 | 6     | L    | 178    |
| 4U2061325007 | 7     | L    | 196    |
| 4U2061326008 | 8     | L    | 214    |
| 4U2061326009 | 9     | L    | 232    |
| 4U2061326010 | 10    | L    | 250    |
| 4U2061326011 | 11    | L    | 268    |
| 4U2061326012 | 12    | L    | 286    |
| 4U2061326013 | 13    | L    | 304    |

R

Screws  
Ø 5.0 mm – Locking Self Tapping Screw  
Ø 4.5 mm – Cortical Self Tapping Screw  
Ø 5.0 mm – Locking Cancellous Screw



R  
Thickness: 5,5  
Width: 16


# FEMUR HIGH OSTEOTOMY PLATE 4.5/5.0 MM



| Code         | Size | Hole | Length |
|--------------|------|------|--------|
| 4U2031915007 | R    | 4    | 113,5  |
| 4U2031925008 | L    | 4    | 113,5  |

Screws  
Ø 5.0 mm – Locking Self Tapping Screw  
Ø 4.5 mm – Cortical Self Tapping Screw  
Ø 5.0 mm – Locking Cancellous Screw

## TIBIA DISTAL MEDIAL LOCKING PLATE



**R**

Thickness: 3,2  
Width: 14

### Screws

Ø 3.5 mm – Locking Self Tapping Screw  
Ø 3.5 mm – Cortical Self Tapping Screw  
Ø 4,0 mm – Locking Cancellous Screw Ø  
4.0 mm – Cancellous Screw Ø 5.0 mm –  
Locking Self Tapping Screw Ø 4.5 mm –  
Cortical Self Tapping Screw Ø 5.0 mm –  
Locking Cancellous Screw

### Code Holes Side Length

|              |    |   |     |
|--------------|----|---|-----|
| 4U2071315003 | 3  | R | 113 |
| 4U2071315004 | 4  | R | 129 |
| 4U2071315005 | 5  | R | 145 |
| 4U2071315006 | 6  | R | 161 |
| 4U2071315007 | 7  | R | 177 |
| 4U2071315008 | 8  | R | 193 |
| 4U2071315009 | 9  | R | 209 |
| 4U2071315010 | 10 | R | 225 |
| 4U2071315011 | 11 | R | 241 |
| 4U2071315012 | 12 | R | 257 |
| 4U2071315013 | 13 | R | 273 |
| 4U2071325003 | 3  | L | 113 |
| 4U2071325004 | 4  | L | 129 |
| 4U2071325005 | 5  | L | 145 |
| 4U2071325006 | 6  | L | 161 |
| 4U2071325007 | 7  | L | 177 |
| 4U2071325008 | 8  | L | 193 |
| 4U2071325009 | 9  | L | 209 |
| 4U2071325010 | 10 | L | 225 |
| 4U2071325011 | 11 | L | 241 |
| 4U2071325012 | 12 | L | 257 |
| 4U2071325013 | 13 | L | 273 |

## TIBIA DISTAL MEDIAL LOCKING PLATE 3.5

### Code Holes Side Length

|              |    |   |     |
|--------------|----|---|-----|
| 4U2071335003 | 3  | R | 117 |
| 4U2071335004 | 4  | R | 130 |
| 4U2071335005 | 5  | R | 143 |
| 4U2071335006 | 6  | R | 156 |
| 4U2071335007 | 7  | R | 169 |
| 4U2071335008 | 8  | R | 182 |
| 4U2071335009 | 9  | R | 195 |
| 4U2071335010 | 10 | R | 208 |
| 4U2071335011 | 11 | R | 221 |
| 4U2071335012 | 12 | R | 234 |
| 4U2071335013 | 13 | R | 247 |
| 4U2071345003 | 3  | L | 117 |
| 4U2071345004 | 4  | L | 130 |
| 4U2071345005 | 5  | L | 143 |
| 4U2071345006 | 6  | L | 156 |
| 4U2071345007 | 7  | L | 169 |
| 4U2071345008 | 8  | L | 182 |
| 4U2071345009 | 9  | L | 195 |
| 4U2071345010 | 10 | L | 208 |
| 4U2071345011 | 11 | L | 221 |
| 4U2071345012 | 12 | L | 234 |
| 4U2071345013 | 13 | L | 247 |

### Screws

Ø 3.5 mm – Locking Self Tapping Screw  
Ø 3.5 mm – Cortical Self Tapping Screw  
Ø 4,0 mm – Locking Cancellous Screw Ø  
4.0 mm – Cancellous Screw



**R**

Thickness: 3,2  
Width: 14

## LARGE NARROW LOCKING PLATE



Thickness: 4,5  
Width: 14

### Code Holes Side Length

|              |    |   |     |
|--------------|----|---|-----|
| 4U2021205004 | 4  | - | 88  |
| 4U2021205005 | 5  | - | 106 |
| 4U2021205006 | 6  | - | 124 |
| 4U2021205007 | 7  | - | 142 |
| 4U2021205008 | 8  | - | 160 |
| 4U2021205009 | 9  | - | 178 |
| 4U2021205010 | 10 | - | 196 |
| 4U2021205011 | 11 | - | 214 |
| 4U2021205012 | 12 | - | 232 |
| 4U2021205013 | 13 | - | 250 |
| 4U2021205014 | 14 | - | 268 |
| 4U2021205015 | 15 | - | 286 |
| 4U2021205016 | 16 | - | 304 |
| 4U2021205018 | 18 | - | 340 |
| 4U2021205020 | 20 | - | 376 |

### Screws

Ø 5.0 mm – Locking Self Tapping Screw  
Ø 4.5 mm – Cortical Self Tapping Screw  
Ø 5.0 mm – Locking Cancellous Screw

## LARGE BROAD LOCKING PLATE



Thickness: 5,5  
Width: 18

### Screws

- Ø 5.0 mm – Locking Self Tapping Screw
- Ø 4.5 mm – Cortical Self Tapping Screw
- Ø 5.0 mm – Locking Cancellous Screw

### Code Holes Side Length

|              |    |   |     |
|--------------|----|---|-----|
| 4U2061205005 | 5  | - | 105 |
| 4U2061205006 | 6  | - | 123 |
| 4U2061205007 | 7  | - | 141 |
| 4U2061205008 | 8  | - | 159 |
| 4U2061205009 | 9  | - | 177 |
| 4U2061205010 | 10 | - | 195 |
| 4U2061205011 | 11 | - | 213 |
| 4U2061205012 | 12 | - | 231 |
| 4U2061205013 | 13 | - | 249 |
| 4U2061205014 | 14 | - | 267 |
| 4U2061205015 | 15 | - | 285 |
| 4U2061205016 | 16 | - | 303 |
| 4U2061205018 | 18 | - | 339 |
| 4U2061205020 | 20 | - | 375 |

## LARGE BROAD LOCKING T PLATE

### Code Holes Side Length

|              |   |   |     |
|--------------|---|---|-----|
| 4U2071105003 | 3 | - | 73  |
| 4U2071105004 | 4 | - | 89  |
| 4U2071105005 | 5 | - | 105 |
| 4U2071105006 | 6 | - | 121 |
| 4U2071105007 | 7 | - | 137 |
| 4U2071105008 | 8 | - | 153 |

### Screws

- Ø 5.0 mm – Locking Self Tapping Screw
- Ø 4.5 mm – Cortical Self Tapping Screw
- Ø 5.0 mm – Locking Cancellous Screw



Thickness: 3  
Width: 15,5

## LARGE BROAD T-SUPPORT LOCKING PLATE

Thickness: 3,5  
Width: 15,5

### Screws

- Ø 5.0 mm – Locking Self Tapping Screw
- Ø 5.0 mm – Locking Cancellous Screw
- Ø 4.5 mm – Cortical Self Tapping Screw

### Code Holes Side Length

|              |   |   |     |
|--------------|---|---|-----|
| 4U2071105014 | 4 | - | 116 |
| 4U2071105015 | 5 | - | 134 |
| 4U2071105016 | 6 | - | 152 |

## TIBIA HIGH OSTEOTOMY PLATE 4.5/5.0 MM

### Code Size Side Length


|              |   |   |     |
|--------------|---|---|-----|
| 4U2031923506 | L | - | 114 |
| 4U2031913505 | S | - | 112 |

- Ø 5.0 mm – Locking Self Tapping Screw
- Ø 4.5 mm – Cortical Self Tapping Screw
- Ø 5.0 mm – Locking Cancellous Screw



Thickness: 3,5  
Width: 16

## LARGE BROAD L-SUPPORT LOCKING PLATE



**R**

Thickness: 3  
Width: 16

Screws

Ø 5.0 mm – Locking Self Tapping Screw

Ø 5.0 mm – Locking Cancellous Screw

Ø 4.5 mm – Cortical Self Tapping Screw

**Code**

**Holes Side Length**

|              |   |   |     |
|--------------|---|---|-----|
| 4U2071115003 | 3 | R | 73  |
| 4U2071115004 | 4 | R | 89  |
| 4U2071115005 | 5 | R | 105 |
| 4U2071115006 | 6 | R | 121 |
| 4U2071115007 | 7 | R | 137 |
| 4U2071115008 | 8 | R | 153 |
| 4U2071125003 | 3 | L | 73  |
| 4U2071125004 | 4 | L | 89  |
| 4U2071125005 | 5 | L | 105 |
| 4U2071125006 | 6 | L | 121 |
| 4U2071125007 | 7 | L | 137 |
| 4U2071125008 | 8 | L | 153 |

## PROXIMAL MEDIAL/ANTEROMEDIAL TIBIA LOCKING PLATE

**Code**

**Holes Side Length**

|              |   |   |     |
|--------------|---|---|-----|
| 4U1071113504 | 4 | R | 108 |
| 4U1071113505 | 5 | R | 125 |
| 4U1071113506 | 6 | R | 142 |
| 4U1071113507 | 7 | R | 159 |
| 4U1071113508 | 8 | R | 176 |
| 4U1071113509 | 9 | R | 193 |
| 4U1071123504 | 4 | L | 108 |
| 4U1071123505 | 5 | L | 125 |
| 4U1071123506 | 6 | L | 142 |
| 4U1071123507 | 7 | L | 159 |
| 4U1071123508 | 8 | L | 176 |
| 4U1071123509 | 9 | L | 193 |

Screws

Ø 5.0 mm – Locking Self Tapping

Screw Ø 4.5 mm – Cortical Self

Tapping Screw Ø 5.0 mm – Locking

Cancellous Screw



**R**

Thickness: 3,5  
Width: 13,5

## TIBIA PROXIMAL LATERAL LOCKING PLATE



**R**

Thickness: 4,5  
Width: 15

Screws

Ø 5.0 mm – Locking Self Tapping Screw

Ø 5.0 mm – Locking Cancellous Screw

Ø 4.5 mm – Cortical Self Tapping Screw


**Code**

**Holes Side Length**

|              |    |   |     |
|--------------|----|---|-----|
| 4U2071113503 | 3  | R | 112 |
| 4U2071113504 | 4  | R | 131 |
| 4U2071113505 | 5  | R | 150 |
| 4U2071113506 | 6  | R | 169 |
| 4U2071113507 | 7  | R | 188 |
| 4U2071113508 | 8  | R | 207 |
| 4U2071113509 | 9  | R | 226 |
| 4U2071113510 | 10 | R | 245 |
| 4U2071113511 | 11 | R | 264 |
| 4U2071113512 | 12 | R | 283 |
| 4U2071113513 | 13 | R | 302 |
| 4U2071123503 | 3  | L | 112 |
| 4U2071123504 | 4  | L | 131 |
| 4U2071123505 | 5  | L | 150 |
| 4U2071123506 | 6  | L | 169 |
| 4U2071123507 | 7  | L | 188 |
| 4U2071123508 | 8  | L | 207 |
| 4U2071123509 | 9  | L | 226 |
| 4U2071125010 | 10 | L | 245 |
| 4U2071125011 | 11 | L | 264 |
| 4U2071125012 | 12 | L | 283 |
| 4U2071125013 | 13 | L | 302 |

## TIBIA DISTAL LATERAL LOCKING PLATE

**R**  
Thickness: 3,8  
Width: 13



### Screws

- Ø 3.5 mm – Locking Self Tapping Screw
- Ø 3.5 mm – Cortical Self Tapping Screw
- Ø 4,0 mm – Locking Cancellous Screw
- Ø 4.0 mm – Cancellous Screw

| Code         | Holes | Side | Length |
|--------------|-------|------|--------|
| 4U1071313503 | 3     | R    | 71     |
| 4U1071313504 | 4     | R    | 84     |
| 4U1071313505 | 5     | R    | 97     |
| 4U1071313506 | 6     | R    | 110    |
| 4U1071313507 | 7     | R    | 123    |
| 4U1071313508 | 8     | R    | 136    |
| 4U1071313509 | 9     | R    | 149    |
| 4U1071313510 | 10    | R    | 162    |
| 4U1071313511 | 11    | R    | 175    |
| 4U1071313512 | 12    | R    | 188    |
| 4U1071313513 | 13    | R    | 201    |
| 4U1071323503 | 3     | L    | 71     |
| 4U1071323504 | 4     | L    | 84     |
| 4U1071323505 | 5     | L    | 97     |
| 4U1071323506 | 6     | L    | 110    |
| 4U1071323507 | 7     | L    | 123    |
| 4U1071323508 | 8     | L    | 136    |
| 4U1071323509 | 9     | L    | 149    |
| 4U1071323510 | 10    | L    | 162    |
| 4U1071323511 | 11    | L    | 175    |
| 4U1071323512 | 12    | L    | 188    |
| 4U1071323513 | 13    | L    | 201    |

## TIBIA DISTAL LATERAL LOCKING PLATE TP2

**Code Holes Side Length**

|              |    |   |     |
|--------------|----|---|-----|
| 4U1071333503 | 3  | R | 71  |
| 4U1071333504 | 4  | R | 84  |
| 4U1071333505 | 5  | R | 97  |
| 4U1071333506 | 6  | R | 110 |
| 4U1071333507 | 7  | R | 123 |
| 4U1071333508 | 8  | R | 136 |
| 4U1071333509 | 9  | R | 149 |
| 4U1071333510 | 10 | R | 162 |
| 4U1071333511 | 11 | R | 175 |
| 4U1071333512 | 12 | R | 188 |
| 4U1071333513 | 13 | R | 201 |
| 4U1071343503 | 3  | L | 71  |
| 4U1071343504 | 4  | L | 84  |
| 4U1071343505 | 5  | L | 97  |
| 4U1071343506 | 6  | L | 110 |
| 4U1071343507 | 7  | L | 123 |
| 4U1071343508 | 8  | L | 136 |
| 4U1071343509 | 9  | L | 149 |
| 4U1071343510 | 10 | L | 162 |
| 4U1071343511 | 11 | L | 175 |
| 4U1071343512 | 12 | L | 188 |
| 4U1071343513 | 13 | L | 201 |

### Screws


- Ø 3.5 mm – Locking Self Tapping Screw
- Ø 3.5 mm – Cortical Self Tapping Screw
- Ø 4,0 mm – Locking Cancellous Screw
- Ø 4.0 mm – Cancellous Screw

**R**  
Thickness: 3,8  
Width: 13



## PILON LOCKING PLATE

Thickness: 3,1  
Width: 11




### Screws

- Ø 3.5 mm – Locking Self Tapping Screw
- Ø 3.5 mm – Cortical Self Tapping Screw
- Ø 4,0 mm – Locking Cancellous Screw
- Ø 4.0 mm – Cancellous Screw

| Code         | Holes | Length |
|--------------|-------|--------|
| 4U1071203515 | 5     | 110    |
| 4U1071203516 | 6     | 125    |
| 4U1071203517 | 7     | 140    |
| 4U1071203518 | 8     | 158    |
| 4U1071203519 | 9     | 170    |

# FIBULA LOCKING PLATE



**R**

Thickness: 2,1  
Width: 9,5

Screws  
Ø 3.5 mm – Locking Self Tapping Screw  
Ø 3.5 mm – Cortical Self Tapping Screw  
Ø 4,0 mm – Locking Cancellous Screw  
Ø 4.0 mm – Cancellous Screw

| Code         | Holes | Side | Length |
|--------------|-------|------|--------|
| 4U1081613502 | 2     | R    | 50     |
| 4U1081613503 | 3     | R    | 62     |
| 4U1081613504 | 4     | R    | 74     |
| 4U1081613505 | 5     | R    | 86     |
| 4U1081613506 | 6     | R    | 98     |
| 4U1081613507 | 7     | R    | 110    |
| 4U1081613508 | 8     | R    | 122    |
| 4U1081613509 | 9     | R    | 134    |
| 4U1081613510 | 10    | R    | 146    |
| 4U1081613511 | 11    | R    | 158    |
| 4U1081613512 | 12    | R    | 170    |
| 4U1081613513 | 13    | R    | 182    |
| 4U1081693502 | 2     | L    | 50     |
| 4U1081693503 | 3     | L    | 62     |
| 4U1081693504 | 4     | L    | 74     |
| 4U1081693505 | 5     | L    | 86     |
| 4U1081693506 | 6     | L    | 98     |
| 4U1081693507 | 7     | L    | 110    |
| 4U1081693508 | 8     | L    | 122    |
| 4U1081693509 | 9     | L    | 134    |
| 4U1081693510 | 10    | L    | 146    |
| 4U1081693511 | 11    | L    | 158    |
| 4U1081693512 | 12    | L    | 170    |
| 4U1081693513 | 13    | L    | 182    |

# FIBULA DISTAL LOCKING PLATE

| Code         | Holes | Side | Length |
|--------------|-------|------|--------|
| 4U1081313502 | 2     | R    | 54     |
| 4U1081313503 | 3     | R    | 66     |
| 4U1081313504 | 4     | R    | 79     |
| 4U1081313505 | 5     | R    | 92     |
| 4U1081313506 | 6     | R    | 105    |
| 4U1081313507 | 7     | R    | 118    |
| 4U1081313508 | 8     | R    | 131    |
| 4U1081313509 | 9     | R    | 144    |
| 4U1081313510 | 10    | R    | 157    |
| 4U1081313511 | 11    | R    | 170    |
| 4U1081313512 | 12    | R    | 183    |
| 4U1081313513 | 13    | R    | 196    |
| 4U1081323502 | 2     | L    | 54     |
| 4U1081323503 | 3     | L    | 66     |
| 4U1081323504 | 4     | L    | 79     |
| 4U1081323505 | 5     | L    | 92     |
| 4U1081323506 | 6     | L    | 105    |
| 4U1081323507 | 7     | L    | 118    |
| 4U1081323508 | 8     | L    | 131    |
| 4U1081323509 | 9     | L    | 144    |
| 4U1081323510 | 10    | L    | 157    |
| 4U1081323511 | 11    | L    | 170    |
| 4U1081323512 | 12    | L    | 183    |
| 4U1081323513 | 13    | L    | 196    |

Screws  
Ø 2.5 mm – Locking Self Tapping Screw  
Ø 2.5 mm – Cortical Self Tapping Screw  
Ø 3.5 mm – Locking Self Tapping Screw  
Ø 3.5 mm – Cortical Self Tapping Screw  
Ø 4,0 mm – Locking Cancellous Screw



**R**

Thickness: 3  
Width: 10

# RECONSTRUCTION SHAFT LOCKING PLATE



Thickness: 3  
Width: 10

Screws  
Ø 3.5 mm – Locking Self Tapping Screw  
Ø 3.5 mm – Cortical Self Tapping Screw  
Ø 4,0 mm – Locking Cancellous Screw

| Code         | Holes | Side | Length |
|--------------|-------|------|--------|
| 4U1071204004 | 4     | -    | 48     |
| 4U1071204005 | 5     | -    | 60     |
| 4U1071204006 | 6     | -    | 72     |
| 4U1071204007 | 7     | -    | 84     |
| 4U1071204008 | 8     | -    | 96     |
| 4U1071204009 | 9     | -    | 108    |
| 4U1071204010 | 10    | -    | 120    |
| 4U1071204011 | 11    | -    | 132    |
| 4U1071204012 | 12    | -    | 144    |
| 4U1071204013 | 13    | -    | 156    |
| 4U1071204014 | 14    | -    | 168    |
| 4U1071204015 | 15    | -    | 180    |
| 4U1071204016 | 16    | -    | 192    |
| 4U1071204017 | 17    | -    | 204    |
| 4U1071204018 | 18    | -    | 216    |
| 4U1071204019 | 19    | -    | 228    |
| 4U1071204020 | 20    | -    | 240    |

## CALCANEUS LOCKING PLATE



| Code         | Holes | Side | Length |
|--------------|-------|------|--------|
| 4U1091013508 | 8     | R    | 60     |
| 4U1091013510 | 10    | R    | 75     |
| 4U1091013515 | 15    | R    | 69     |
| 4U1091013516 | 15    | R    | 75     |
| 4U1091023508 | 8     | L    | 60     |
| 4U1091023510 | 10    | L    | 75     |
| 4U1091023515 | 15    | L    | 69     |
| 4U1091023516 | 15    | L    | 75     |

## SMALL FEMUR PROXIMAL HIP LOCKING PLATE

3.5

| Code         | Holes | Angle | Length |
|--------------|-------|-------|--------|
| 4U1061103503 | 3     | 100   | 76     |
| 4U1061103513 | 3     | 110   | 76     |
| 4U1061103504 | 3     | 120   | 62     |
| 4U1061103506 | 6     | 120   | 114    |
| 4U1061103523 | 3     | 150   | 65     |

5.0

| Code         | Hol | Angle | Length |
|--------------|-----|-------|--------|
| 4U2061105003 | 4   | 100   | 116    |
| 4U2061105013 | 4   | 110   | 117    |
| 4U2061105004 | 4   | 120   | 117    |
| 4U2061105006 | 6   | 120   | 149    |
| 4U2061105023 | 4   | 150   | 117    |

Ø 3.5 mm – Locking Self Tapping Screw  
 Ø 3.5 mm – Cortical Self Tapping Screw  
 Ø 5,0 mm – Locking Cancellous Screw  
 Ø 5,0 mm – Cancellous Screw  
 Ø 4.5 mm – Cortical Self Tapping Screw



## PELVIC RECONSTRUCTION CURVED LOCKING PLATE



| Code         | Holes | Side | Length |
|--------------|-------|------|--------|
| 4U2061003504 | 4     | -    | 47     |
| 4U2061003505 | 5     | -    | 61     |
| 4U2061003506 | 6     | -    | 73     |
| 4U2061003507 | 7     | -    | 85     |
| 4U2061003508 | 8     | -    | 97     |
| 4U2061003509 | 9     | -    | 108    |
| 4U2061003510 | 10    | -    | 119    |
| 4U2061003511 | 11    | -    | 130    |
| 4U2061003512 | 12    | -    | 141    |
| 4U2061003513 | 13    | -    | 150    |
| 4U2061003514 | 14    | -    | 157    |

Ø 3.5 mm – Locking Self Tapping Screw  
 Ø 3.5 mm – Cortical Self Tapping Screw  
 Ø 4,0 mm – Locking Cancellous Screw

## PELVIC SYMPHYSEAL PLATE



Thickness: 2,5  
Width: 10

### Screws

Ø 3.5 mm – Locking Self Tapping Screw

Ø 3.5 mm – Cortical Self Tapping Screw

Ø 4,0 mm – Locking Cancellous Screw

| Code         | Holes | Side | Length |
|--------------|-------|------|--------|
| 4U2061003552 | 4     | -    | 52     |
| 4U2061003555 | 4     | -    | 55     |
| 4U2061003558 | 4     | -    | 58     |

## PELVIC RECONSTRUCTION STRAIGHT LOCKING PLATE

| Code         | Holes | Side | Length |
|--------------|-------|------|--------|
| 4U2061013504 | 4     | -    | 46     |
| 4U2061013505 | 5     | -    | 58     |
| 4U2061013506 | 6     | -    | 70     |
| 4U2061013507 | 7     | -    | 82     |
| 4U2061013508 | 8     | -    | 94     |
| 4U2061013509 | 9     | -    | 106    |
| 4U2061013510 | 10    | -    | 118    |
| 4U2061013511 | 11    | -    | 130    |
| 4U2061013512 | 12    | -    | 142    |
| 4U2061013513 | 13    | -    | 154    |
| 4U2061013514 | 14    | -    | 166    |
| 4U2061013515 | 15    | -    | 178    |
| 4U2061013516 | 16    | -    | 190    |



Thickness: 2,5  
Width: 10

### Screws

Ø 3.5 mm – Locking Self Tapping Screw

Ø 3.5 mm – Cortical Self Tapping Screw

Ø 4,0 mm – Locking Cancellous Screw

## EIGHT NON LOCKING PLATE



Thickness: 1,6  
Width: 9,8

### Screws

Ø 3.5 mm – Cortical Self Tapping Screw

Ø 3.5 mm – Cannulated Screw

Ø 4,0 mm – Cancellous Screw

Ø 4,5 mm – Full Tread Cannulated Screw

| Code         | Size | Side | Length |
|--------------|------|------|--------|
| 4U8880000050 | STD. | -    | 26     |
| 4U8880000051 | LONG | -    | 30     |

**LOCKING SELF TAPPING  
CORTICAL SCREW Ø 1,5 MM**



| Code         | Size  |
|--------------|-------|
| 4U3001001506 | 6 mm  |
| 4U3001001507 | 7 mm  |
| 4U3001001508 | 8 mm  |
| 4U3001001509 | 9 mm  |
| 4U3001001510 | 10 mm |
| 4U3001001511 | 11 mm |
| 4U3001001512 | 12 mm |
| 4U3001001513 | 13 mm |
| 4U3001001514 | 14 mm |
| 4U3001001515 | 15 mm |
| 4U3001001516 | 16 mm |
| 4U3001001517 | 17 mm |
| 4U3001001518 | 18 mm |
| 4U3001001519 | 19 mm |
| 4U3001001520 | 20 mm |
| 4U3001001522 | 22 mm |
| 4U3001001524 | 24 mm |
| 4U3001001526 | 26 mm |
| 4U3001001528 | 28 mm |
| 4U3001001530 | 30 mm |

**LOCKING SELF TAPPING  
CORTICAL SCREW Ø 2,0 MM**



| Code         | Size  |
|--------------|-------|
| 4U3001002006 | 6 mm  |
| 4U3001002007 | 7 mm  |
| 4U3001002008 | 8 mm  |
| 4U3001002009 | 9 mm  |
| 4U3001002010 | 10 mm |
| 4U3001002011 | 11 mm |
| 4U3001002012 | 12 mm |
| 4U3001002013 | 13 mm |
| 4U3001002014 | 14 mm |
| 4U3001002015 | 15 mm |
| 4U3001002016 | 16 mm |
| 4U3001002017 | 17 mm |
| 4U3001002018 | 18 mm |
| 4U3001002019 | 19 mm |
| 4U3001002020 | 20 mm |
| 4U3001002022 | 22 mm |
| 4U3001002024 | 24 mm |
| 4U3001002026 | 26 mm |
| 4U3001002028 | 28 mm |
| 4U3001002030 | 30 mm |

**LOCKING SELF TAPPING  
CORTICAL SCREW Ø 2,5 MM**



| Code         | Size  |
|--------------|-------|
| 4U3001002506 | 6 mm  |
| 4U3001002507 | 7 mm  |
| 4U3001002508 | 8 mm  |
| 4U3001002509 | 9 mm  |
| 4U3001002510 | 10 mm |
| 4U3001002511 | 11 mm |
| 4U3001002512 | 12 mm |
| 4U3001002513 | 13 mm |
| 4U3001002514 | 14 mm |
| 4U3001002515 | 15 mm |
| 4U3001002516 | 16 mm |
| 4U3001002517 | 17 mm |
| 4U3001002518 | 18 mm |
| 4U3001002519 | 19 mm |
| 4U3001002520 | 20 mm |
| 4U3001002522 | 22 mm |
| 4U3001002524 | 24 mm |
| 4U3001002526 | 26 mm |
| 4U3001002528 | 28 mm |
| 4U3001002530 | 30 mm |
| 4U3001002532 | 32 mm |
| 4U3001002534 | 34 mm |
| 4U3001002536 | 36 mm |
| 4U3001002538 | 38 mm |
| 4U3001002540 | 40 mm |

**LOCKING SELF TAPPING  
CORTICAL SCREW Ø 2,7 MM**



**LOCKING SELF TAPPING  
CORTICAL SCREW Ø 3,5 MM**



**LOCKING SELF TAPPING  
CORTICAL SCREW Ø 5,0 MM**



| Code         | Size  | Code         | Size  | Code         | Size  | Code         | Size   |
|--------------|-------|--------------|-------|--------------|-------|--------------|--------|
| 4U3001002710 | 10 mm | 4U3001003510 | 10 mm | 4U3001005016 | 16 mm | 4U3001005065 | 65 mm  |
| 4U3001002711 | 11 mm | 4U3001003512 | 12 mm | 4U3001005018 | 18 mm | 4U3001005066 | 66 mm  |
| 4U3001002712 | 12 mm | 4U3001003514 | 14 mm | 4U3001005020 | 20 mm | 4U3001005068 | 68 mm  |
| 4U3001002713 | 13 mm | 4U3001003516 | 16 mm | 4U3001005022 | 22 mm | 4U3001005070 | 70 mm  |
| 4U3001002714 | 14 mm | 4U3001003518 | 18 mm | 4U3001005024 | 24 mm | 4U3001005072 | 72 mm  |
| 4U3001002715 | 15 mm | 4U3001003520 | 20 mm | 4U3001005026 | 26 mm | 4U3001005074 | 74 mm  |
| 4U3001002716 | 16 mm | 4U3001003522 | 22 mm | 4U3001005028 | 28 mm | 4U3001005075 | 75 mm  |
| 4U3001002717 | 17 mm | 4U3001003524 | 24 mm | 4U3001005030 | 30 mm | 4U3001005076 | 76 mm  |
| 4U3001002718 | 18 mm | 4U3001003526 | 26 mm | 4U3001005032 | 32 mm | 4U3001005078 | 78 mm  |
| 4U3001002719 | 19 mm | 4U3001003528 | 28 mm | 4U3001005034 | 34 mm | 4U3001005080 | 80 mm  |
| 4U3001002720 | 20 mm | 4U3001003530 | 30 mm | 4U3001005036 | 36 mm | 4U3001005082 | 82 mm  |
| 4U3001002722 | 22 mm | 4U3001003532 | 32 mm | 4U3001005038 | 38 mm | 4U3001005084 | 84 mm  |
| 4U3001002724 | 24 mm | 4U3001003534 | 34 mm | 4U3001005040 | 40 mm | 4U3001005085 | 85 mm  |
| 4U3001002726 | 26 mm | 4U3001003536 | 36 mm | 4U3001005042 | 42 mm | 4U3001005086 | 86 mm  |
| 4U3001002728 | 28 mm | 4U3001003538 | 38 mm | 4U3001005044 | 44 mm | 4U3001005088 | 88 mm  |
| 4U3001002730 | 30 mm | 4U3001003540 | 40 mm | 4U3001005046 | 46 mm | 4U3001005090 | 90 mm  |
| 4U3001002732 | 32 mm | 4U3001003542 | 42 mm | 4U3001005048 | 48 mm | 4U3001005092 | 92 mm  |
| 4U3001002734 | 34 mm | 4U3001003544 | 44 mm | 4U3001005050 | 50 mm | 4U3001005094 | 94 mm  |
| 4U3001002736 | 36 mm | 4U3001003546 | 46 mm | 4U3001005052 | 52 mm | 4U3001005096 | 96 mm  |
| 4U3001002738 | 38 mm | 4U3001003548 | 48 mm | 4U3001005054 | 54 mm | 4U3001005098 | 98 mm  |
| 4U3001002740 | 40 mm | 4U3001003550 | 50 mm | 4U3001005055 | 55 mm | 4U3001005010 | 100 mm |
| 4U3001002742 | 42 mm | 4U3001003552 | 52 mm | 4U3001005056 | 56 mm | 4U3001005011 | 102 mm |
| 4U3001002744 | 44 mm | 4U3001003554 | 54 mm | 4U3001005058 | 58 mm | 4U3001005012 | 104 mm |
| 4U3001002746 | 46 mm | 4U3001003555 | 55 mm | 4U3001005060 | 60 mm | 4U3001005013 | 106 mm |
| 4U3001002748 | 48 mm | 4U3001003556 | 56 mm | 4U3001005062 | 62 mm | 4U3001005014 | 108 mm |
| 4U3001002750 | 50 mm | 4U3001003558 | 58 mm | 4U3001005064 | 64 mm | 4U3001005015 | 110 mm |
|              |       | 4U3001003560 | 60 mm |              |       |              |        |
|              |       | 4U3001003565 | 65 mm |              |       |              |        |
|              |       | 4U3001003570 | 70 mm |              |       |              |        |

**CORTICAL SELF TAPPING SCREW**  
**Ø 1,5 MM**



**CORTICAL SELF TAPPING SCREW**  
**Ø 2,0 MM**



**CORTICAL SELF TAPPING SCREW**  
**Ø 2,5 MM**



**Code Size**

|              |       |
|--------------|-------|
| 4U3002001506 | 6 mm  |
| 4U3002001507 | 7 mm  |
| 4U3002001508 | 8 mm  |
| 4U3002001509 | 9 mm  |
| 4U3002001510 | 10 mm |
| 4U3002001511 | 11 mm |
| 4U3002001512 | 12 mm |
| 4U3002001513 | 13 mm |
| 4U3002001514 | 14 mm |
| 4U3002001515 | 15 mm |
| 4U3002001516 | 16 mm |
| 4U3002001517 | 17 mm |
| 4U3002001518 | 18 mm |
| 4U3002001519 | 19 mm |
| 4U3002001520 | 20 mm |
| 4U3002001522 | 22 mm |
| 4U3002001524 | 24 mm |
| 4U3002001526 | 26 mm |
| 4U3002001528 | 28 mm |
| 4U3002001530 | 30 mm |

**Code Size**

|              |       |
|--------------|-------|
| 4U3002002006 | 6 mm  |
| 4U3002002007 | 7 mm  |
| 4U3002002008 | 8 mm  |
| 4U3002002009 | 9 mm  |
| 4U3002002010 | 10 mm |
| 4U3002002011 | 11 mm |
| 4U3002002012 | 12 mm |
| 4U3002002013 | 13 mm |
| 4U3002002014 | 14 mm |
| 4U3002002015 | 15 mm |
| 4U3002002016 | 16 mm |
| 4U3002002017 | 17 mm |
| 4U3002002018 | 18 mm |
| 4U3002002019 | 19 mm |
| 4U3002002020 | 20 mm |
| 4U3002002022 | 22 mm |
| 4U3002002024 | 24 mm |
| 4U3002002026 | 26 mm |
| 4U3002002028 | 28 mm |
| 4U3002002030 | 30 mm |

**Code Size**

|              |       |
|--------------|-------|
| 4U3002002506 | 6 mm  |
| 4U3002002507 | 7 mm  |
| 4U3002002508 | 8 mm  |
| 4U3002002509 | 9 mm  |
| 4U3002002510 | 10 mm |
| 4U3002002511 | 11 mm |
| 4U3002002512 | 12 mm |
| 4U3002002513 | 13 mm |
| 4U3002002514 | 14 mm |
| 4U3002002515 | 15 mm |
| 4U3002002516 | 16 mm |
| 4U3002002517 | 17 mm |
| 4U3002002518 | 18 mm |
| 4U3002002519 | 19 mm |
| 4U3002002520 | 20 mm |
| 4U3002002522 | 22 mm |
| 4U3002002524 | 24 mm |
| 4U3002002526 | 26 mm |
| 4U3002002528 | 28 mm |
| 4U3002002530 | 30 mm |
| 4U3002002532 | 32 mm |
| 4U3002002534 | 34 mm |
| 4U3002002536 | 36 mm |
| 4U3002002538 | 38 mm |
| 4U3002002540 | 40 mm |

**CORTICAL SELF TAPPING SCREW**  
**Ø 2,7 MM**

| Code         | Size  |
|--------------|-------|
| 4U3002002710 | 10 mm |
| 4U3002002712 | 12 mm |
| 4U3002002714 | 14 mm |
| 4U3002002716 | 16 mm |
| 4U3002002718 | 18 mm |
| 4U3002002720 | 20 mm |
| 4U3002002722 | 22 mm |
| 4U3002002724 | 24 mm |
| 4U3002002726 | 26 mm |
| 4U3002002728 | 28 mm |
| 4U3002002730 | 30 mm |
| 4U3002002732 | 32 mm |
| 4U3002002734 | 34 mm |
| 4U3002002736 | 36 mm |
| 4U3002002738 | 38 mm |
| 4U3002002740 | 40 mm |
| 4U3002002742 | 42 mm |
| 4U3002002744 | 44 mm |
| 4U3002002746 | 46 mm |
| 4U3002002748 | 48 mm |
| 4U3002002750 | 50 mm |

**CORTICAL SELF TAPPING SCREW**  
**Ø 3,5 MM**

| Code         | Size  |
|--------------|-------|
| 4U3002003510 | 10 mm |
| 4U3002003512 | 12 mm |
| 4U3002003514 | 14 mm |
| 4U3002003516 | 16 mm |
| 4U3002003518 | 18 mm |
| 4U3002003520 | 20 mm |
| 4U3002003522 | 22 mm |
| 4U3002003524 | 24 mm |
| 4U3002003526 | 26 mm |
| 4U3002003528 | 28 mm |
| 4U3002003530 | 30 mm |
| 4U3002003532 | 32 mm |
| 4U3002003534 | 34 mm |
| 4U3002003536 | 36 mm |
| 4U3002003538 | 38 mm |
| 4U3002003540 | 40 mm |
| 4U3002003542 | 42 mm |
| 4U3002003544 | 44 mm |
| 4U3002003545 | 45 mm |
| 4U3002003546 | 46 mm |
| 4U3002003548 | 48 mm |
| 4U3002003550 | 50 mm |
| 4U3002003555 | 55 mm |
| 4U3002003560 | 60 mm |
| 4U3002003565 | 65 mm |
| 4U3002003570 | 70 mm |
| 4U3002003575 | 75 mm |
| 4U3002003580 | 80 mm |
| 4U3002003585 | 85 mm |
| 4U3002003590 | 90 mm |

**CORTICAL SELF TAPPING SCREW**  
**Ø 4,5 MM**

| Code         | Size  |
|--------------|-------|
| 4U3002004514 | 14 mm |
| 4U3002004516 | 16 mm |
| 4U3002004518 | 18 mm |
| 4U3002004520 | 20 mm |
| 4U3002004522 | 22 mm |
| 4U3002004524 | 24 mm |
| 4U3002004526 | 26 mm |
| 4U3002004528 | 28 mm |
| 4U3002004530 | 30 mm |
| 4U3002004532 | 32 mm |
| 4U3002004534 | 34 mm |
| 4U3002004536 | 36 mm |
| 4U3002004538 | 38 mm |
| 4U3002004540 | 40 mm |
| 4U3002004542 | 42 mm |
| 4U3002004544 | 44 mm |
| 4U3002004546 | 46 mm |
| 4U3002004548 | 48 mm |
| 4U3002004550 | 50 mm |
| 4U3002004552 | 52 mm |
| 4U3002004554 | 54 mm |
| 4U3002004555 | 55 mm |
| 4U3002004556 | 56 mm |
| 4U3002004558 | 58 mm |
| 4U3002004560 | 60 mm |
| 4U3002004562 | 62 mm |
| 4U3002004564 | 64 mm |
| 4U3002004565 | 65 mm |
| 4U3002004566 | 66 mm |
| 4U3002004568 | 68 mm |
| 4U3002004570 | 70 mm |

**LOCKING CANCELLOUS/  
SPONGIOSE SCREW Ø 4,0 MM**



| Code         | Size  |
|--------------|-------|
| 4U3008034010 | 10 mm |
| 4U3008034012 | 12 mm |
| 4U3008034014 | 14 mm |
| 4U3008034016 | 16 mm |
| 4U3008034018 | 18 mm |
| 4U3008034020 | 20 mm |
| 4U3008034022 | 22 mm |
| 4U3008034024 | 24 mm |
| 4U3008034026 | 26 mm |
| 4U3008034028 | 28 mm |
| 4U3008034030 | 30 mm |
| 4U3008034035 | 35 mm |
| 4U3008034040 | 40 mm |
| 4U3008034045 | 45 mm |
| 4U3008034050 | 50 mm |
| 4U3008034055 | 55 mm |
| 4U3008034060 | 60 mm |
| 4U3008034065 | 65 mm |
| 4U3008034070 | 70 mm |
| 4U3008034075 | 75 mm |

**LOCKING CANCELLOUS/  
SPONGIOSE SCREW Ø 5,0 MM**



| Code         | Size  | Code         | Size   |
|--------------|-------|--------------|--------|
| 4U3008035020 | 20 mm | 4U3008035066 | 66 mm  |
| 4U3008035022 | 22 mm | 4U3008035068 | 68 mm  |
| 4U3008035024 | 24 mm | 4U3008035070 | 70 mm  |
| 4U3008035026 | 26 mm | 4U3008035072 | 72 mm  |
| 4U3008035028 | 28 mm | 4U3008035074 | 74 mm  |
| 4U3008035030 | 30 mm | 4U3008035075 | 75 mm  |
| 4U3008035032 | 32 mm | 4U3008035076 | 76 mm  |
| 4U3008035034 | 34 mm | 4U3008035078 | 78 mm  |
| 4U3008035036 | 36 mm | 4U3008035080 | 80 mm  |
| 4U3008035038 | 38 mm | 4U3008035082 | 82 mm  |
| 4U3008035040 | 40 mm | 4U3008035084 | 84 mm  |
| 4U3008035042 | 42 mm | 4U3008035085 | 85 mm  |
| 4U3008035044 | 44 mm | 4U3008035086 | 86 mm  |
| 4U3008035046 | 46 mm | 4U3008035088 | 88 mm  |
| 4U3008035048 | 48 mm | 4U3008035090 | 90 mm  |
| 4U3008035050 | 50 mm | 4U3008035092 | 92 mm  |
| 4U3008035052 | 52 mm | 4U3008035094 | 94 mm  |
| 4U3008035054 | 54 mm | 4U3008035096 | 96 mm  |
| 4U3008035055 | 55 mm | 4U3008035098 | 98 mm  |
| 4U3008035056 | 56 mm | 4U3008035010 | 100 mm |
| 4U3008035058 | 58 mm | 4U3008035011 | 102 mm |
| 4U3008035060 | 60 mm | 4U3008035012 | 104 mm |
| 4U3008035062 | 62 mm | 4U3008035013 | 106 mm |
| 4U3008035064 | 64 mm | 4U3008035014 | 108 mm |
| 4U3008035065 | 65 mm | 4U3008035015 | 110 mm |

**CANCELLOUS/SPONGIOSE SCREW FULL THREAD  
Ø 4,0 MM**



| Code         | Size  |
|--------------|-------|
| 4U3003034010 | 10 mm |
| 4U3003034012 | 12 mm |
| 4U3003034014 | 14 mm |
| 4U3003034016 | 16 mm |
| 4U3003034018 | 18 mm |
| 4U3003034020 | 20 mm |
| 4U3003034022 | 22 mm |
| 4U3003034024 | 24 mm |
| 4U3003034026 | 26 mm |
| 4U3003034028 | 28 mm |
| 4U3003034030 | 30 mm |
| 4U3003034032 | 32 mm |
| 4U3003034034 | 34 mm |
| 4U3003034036 | 36 mm |
| 4U3003034038 | 38 mm |
| 4U3003034040 | 40 mm |
| 4U3003034045 | 45 mm |
| 4U3003034050 | 50 mm |
| 4U3003034055 | 55 mm |
| 4U3003034060 | 60 mm |
| 4U3003034065 | 65 mm |
| 4U3003034070 | 70 mm |

**CANCELLOUS/SPONGIOSE SCREW FULL THREAD  
Ø 6,5 MM**



| Code         | Size   |
|--------------|--------|
| 4U3003036530 | 30 mm  |
| 4U3003036535 | 35 mm  |
| 4U3003036540 | 40 mm  |
| 4U3003036545 | 45 mm  |
| 4U3003036550 | 50 mm  |
| 4U3003036555 | 55 mm  |
| 4U3003036560 | 60 mm  |
| 4U3003036565 | 65 mm  |
| 4U3003036570 | 70 mm  |
| 4U3003036575 | 75 mm  |
| 4U3003036580 | 80 mm  |
| 4U3003036585 | 85 mm  |
| 4U3003036590 | 90 mm  |
| 4U3003036595 | 95 mm  |
| 4U3003036510 | 100 mm |
| 4U3003036515 | 105 mm |
| 4U3003036511 | 110 mm |
| 4U3003036516 | 115 mm |
| 4U3003036512 | 120 mm |

**CANCELLOUS SCREW HALF  
THREAD Ø 4,0 MM**



| Code         | Size  |
|--------------|-------|
| 4U3003044030 | 30 mm |
| 4U3003044035 | 35 mm |
| 4U3003044040 | 40 mm |
| 4U3003044045 | 45 mm |
| 4U3003044050 | 50 mm |
| 4U3003044055 | 55 mm |
| 4U3003044060 | 60 mm |
| 4U3003044065 | 65 mm |
| 4U3003044070 | 70 mm |

**CANCELLOUS SCREW 16 MM  
HALF THREAD Ø 6,5 MM**



| Code         | Size   |
|--------------|--------|
| 4U3003046545 | 45 mm  |
| 4U3003046550 | 50 mm  |
| 4U3003046555 | 55 mm  |
| 4U3003046560 | 60 mm  |
| 4U3003046565 | 65 mm  |
| 4U3003046570 | 70 mm  |
| 4U3003046575 | 75 mm  |
| 4U3003046580 | 80 mm  |
| 4U3003046585 | 85 mm  |
| 4U3003046590 | 90 mm  |
| 4U3003046595 | 95 mm  |
| 4U3003046510 | 100 mm |
| 4U3003046515 | 105 mm |
| 4U3003046511 | 110 mm |
| 4U3003046516 | 115 mm |
| 4U3003046512 | 120 mm |

**CANCELLOUS SCREW 32 MM  
THREAD Ø 6,5 MM**



| Code         | Size   |
|--------------|--------|
| 4U3003047545 | 45 mm  |
| 4U3003047550 | 50 mm  |
| 4U3003047555 | 55 mm  |
| 4U3003047560 | 60 mm  |
| 4U3003047565 | 65 mm  |
| 4U3003047570 | 70 mm  |
| 4U3003047575 | 75 mm  |
| 4U3003047580 | 80 mm  |
| 4U3003047585 | 85 mm  |
| 4U3003047590 | 90 mm  |
| 4U3003047595 | 95 mm  |
| 4U3003047600 | 100 mm |
| 4U3003047605 | 105 mm |
| 4U3003047610 | 110 mm |
| 4U3003047615 | 115 mm |
| 4U3003047620 | 120 mm |

**CANNULATED COMPRESSION  
SCREW Ø 2,7 MM**



| Code         | Size  |
|--------------|-------|
| 4U3004032720 | 20 mm |
| 4U3004032722 | 22 mm |
| 4U3004032724 | 24 mm |
| 4U3004032726 | 26 mm |
| 4U3004032728 | 28 mm |
| 4U3004032730 | 30 mm |
| 4U3004032732 | 32 mm |
| 4U3004032734 | 34 mm |
| 4U3004032736 | 36 mm |
| 4U3004032738 | 38 mm |
| 4U3004032740 | 40 mm |
| 4U3004032742 | 42 mm |
| 4U3004032744 | 44 mm |
| 4U3004032746 | 46 mm |
| 4U3004032748 | 48 mm |
| 4U3004032750 | 50 mm |

**CANNULATED COMPRESSION  
SCREW FULL THREAD Ø 2,7 MM**



| Code         | Size  |
|--------------|-------|
| 4U3004044520 | 20 mm |
| 4U3004044522 | 22 mm |
| 4U3004044524 | 24 mm |
| 4U3004044526 | 26 mm |
| 4U3004044528 | 28 mm |
| 4U3004044530 | 30 mm |
| 4U3004044532 | 32 mm |
| 4U3004044534 | 34 mm |
| 4U3004044536 | 36 mm |
| 4U3004044538 | 38 mm |
| 4U3004044540 | 40 mm |
| 4U3004044545 | 45 mm |
| 4U3004044550 | 50 mm |
| 4U3004044555 | 55 mm |
| 4U3004044560 | 60 mm |
| 4U3004044565 | 65 mm |
| 4U3004044570 | 70 mm |
| 4U3004044575 | 75 mm |
| 4U3004044580 | 80 mm |

**CANNULATED COMPRESSION  
SCREW Ø 3,5 MM**



| Code         | Size  |
|--------------|-------|
| 4U3004033520 | 20 mm |
| 4U3004033522 | 22 mm |
| 4U3004033524 | 24 mm |
| 4U3004033526 | 26 mm |
| 4U3004033528 | 28 mm |
| 4U3004033530 | 30 mm |
| 4U3004033532 | 32 mm |
| 4U3004033534 | 34 mm |
| 4U3004033536 | 36 mm |
| 4U3004033538 | 38 mm |
| 4U3004033540 | 40 mm |
| 4U3004033542 | 42 mm |
| 4U3004033544 | 44 mm |
| 4U3004033546 | 46 mm |
| 4U3004033548 | 48 mm |
| 4U3004033550 | 50 mm |

**CANNULATED COMPRESSION SCREW**  
Ø 4,5 MM



| Code         | Size  |
|--------------|-------|
| 4U3004034525 | 25 mm |
| 4U3004034530 | 30 mm |
| 4U3004034535 | 35 mm |
| 4U3004034540 | 40 mm |
| 4U3004034545 | 45 mm |
| 4U3004034550 | 50 mm |
| 4U3004034555 | 55 mm |
| 4U3004034560 | 60 mm |
| 4U3004034565 | 65 mm |
| 4U3004034570 | 70 mm |
| 4U3004034575 | 75 mm |
| 4U3004034580 | 80 mm |
| 4U3004034585 | 85 mm |
| 4U3004034590 | 90 mm |

**CANNULATED COMPRESSION SCREW**  
16 MM THREAD Ø 6,5 MM




| Code         | Size   |
|--------------|--------|
| 4U3004036550 | 50 mm  |
| 4U3004036555 | 55 mm  |
| 4U3004036560 | 60 mm  |
| 4U3004036565 | 65 mm  |
| 4U3004036570 | 70 mm  |
| 4U3004036575 | 75 mm  |
| 4U3004036580 | 80 mm  |
| 4U3004036585 | 85 mm  |
| 4U3004036590 | 90 mm  |
| 4U3004036595 | 95 mm  |
| 4U3004036100 | 100 mm |
| 4U3004036105 | 105 mm |
| 4U3004036110 | 110 mm |
| 4U3004036115 | 115 mm |
| 4U3004036120 | 120 mm |

**CANNULATED COMPRESSION SCREW**  
32 MM THREAD Ø 6,5 MM



| Code         | Size   |
|--------------|--------|
| 4U3004036645 | 45 mm  |
| 4U3004036650 | 50 mm  |
| 4U3004036655 | 55 mm  |
| 4U3004036660 | 60 mm  |
| 4U3004036665 | 65 mm  |
| 4U3004036670 | 70 mm  |
| 4U3004036675 | 75 mm  |
| 4U3004036680 | 80 mm  |
| 4U3004036685 | 85 mm  |
| 4U3004036690 | 90 mm  |
| 4U3004036695 | 95 mm  |
| 4U3004036700 | 100 mm |
| 4U3004036705 | 105 mm |
| 4U3004036710 | 110 mm |
| 4U3004036715 | 115 mm |
| 4U3004036720 | 120 mm |

**CANNULATED COMPRESSION SCREW**  
16 MM THREAD Ø 7,3 MM



| Code         | Size   |
|--------------|--------|
| 4U3004037350 | 50 mm  |
| 4U3004037355 | 55 mm  |
| 4U3004037360 | 60 mm  |
| 4U3004037365 | 65 mm  |
| 4U3004037370 | 70 mm  |
| 4U3004037375 | 75 mm  |
| 4U3004037380 | 80 mm  |
| 4U3004037385 | 85 mm  |
| 4U3004037390 | 90 mm  |
| 4U3004037395 | 95 mm  |
| 4U3004037100 | 100 mm |
| 4U3004037105 | 105 mm |
| 4U3004037110 | 110 mm |
| 4U3004037115 | 115 mm |
| 4U3004037120 | 120 mm |

**CANNULATED COMPRESSION SCREW**  
32 MM THREAD Ø 7,3 MM



| Code         | Size   |
|--------------|--------|
| 4U3004047345 | 45 mm  |
| 4U3004047350 | 50 mm  |
| 4U3004047355 | 55 mm  |
| 4U3004047360 | 60 mm  |
| 4U3004047365 | 65 mm  |
| 4U3004047370 | 70 mm  |
| 4U3004047375 | 75 mm  |
| 4U3004047380 | 80 mm  |
| 4U3004047385 | 85 mm  |
| 4U3004047390 | 90 mm  |
| 4U3004047395 | 95 mm  |
| 4U3004047400 | 100 mm |
| 4U3004047405 | 105 mm |
| 4U3004047410 | 110 mm |
| 4U3004047415 | 115 mm |

**LOCKING CANNULATED  
FULL THREAD SCREW Ø 2,5 MM**



**Code Size**

|              |       |
|--------------|-------|
| 4U3005032506 | 6 mm  |
| 4U3005032507 | 7 mm  |
| 4U3005032508 | 8 mm  |
| 4U3005032509 | 9 mm  |
| 4U3005032510 | 10 mm |
| 4U3005032511 | 11 mm |
| 4U3005032512 | 12 mm |
| 4U3005032513 | 13 mm |
| 4U3005032514 | 14 mm |
| 4U3005032515 | 15 mm |
| 4U3005032516 | 16 mm |
| 4U3005032517 | 17 mm |
| 4U3005032518 | 18 mm |
| 4U3005032519 | 19 mm |
| 4U3005032520 | 20 mm |
| 4U3005032522 | 22 mm |
| 4U3005032524 | 24 mm |
| 4U3005032526 | 26 mm |
| 4U3005032528 | 28 mm |
| 4U3005032530 | 30 mm |

**LOCKING CANNULATED  
FULL THREAD SCREW Ø 3,5 MM**



**Code Size**

|              |       |
|--------------|-------|
| 4U3005033510 | 10 mm |
| 4U3005033512 | 12 mm |
| 4U3005033514 | 14 mm |
| 4U3005033516 | 16 mm |
| 4U3005033518 | 18 mm |
| 4U3005033520 | 20 mm |
| 4U3005033522 | 22 mm |
| 4U3005033524 | 24 mm |
| 4U3005033526 | 26 mm |
| 4U3005033528 | 28 mm |
| 4U3005033530 | 30 mm |
| 4U3005033532 | 32 mm |
| 4U3005033534 | 34 mm |
| 4U3005033536 | 36 mm |
| 4U3005033538 | 38 mm |
| 4U3005033540 | 40 mm |
| 4U3005033542 | 42 mm |
| 4U3005033544 | 44 mm |
| 4U3005033546 | 46 mm |
| 4U3005033548 | 48 mm |
| 4U3005033550 | 50 mm |
| 4U3005033554 | 54 mm |
| 4U3005033555 | 55 mm |
| 4U3005033556 | 56 mm |
| 4U3005033558 | 58 mm |
| 4U3005033560 | 60 mm |

**LOCKING CANNULATED FULL  
THREAD SCREW Ø 5,0 MM**



**Code Size Code Size**

|              |       |              |        |
|--------------|-------|--------------|--------|
| 4U3005035016 | 16 mm | 4U3005035065 | 65 mm  |
| 4U3005035018 | 18 mm | 4U3005035066 | 66 mm  |
| 4U3005035020 | 20 mm | 4U3005035068 | 68 mm  |
| 4U3005035022 | 22 mm | 4U3005035070 | 70 mm  |
| 4U3005035024 | 24 mm | 4U3005035072 | 72 mm  |
| 4U3005035026 | 26 mm | 4U3005035074 | 74 mm  |
| 4U3005035028 | 28 mm | 4U3005035075 | 75 mm  |
| 4U3005035030 | 30 mm | 4U3005035076 | 76 mm  |
| 4U3005035032 | 32 mm | 4U3005035078 | 78 mm  |
| 4U3005035034 | 34 mm | 4U3005035080 | 80 mm  |
| 4U3005035036 | 36 mm | 4U3005035082 | 82 mm  |
| 4U3005035038 | 38 mm | 4U3005035084 | 84 mm  |
| 4U3005035040 | 40 mm | 4U3005035085 | 85 mm  |
| 4U3005035042 | 42 mm | 4U3005035086 | 86 mm  |
| 4U3005035044 | 44 mm | 4U3005035088 | 88 mm  |
| 4U3005035046 | 46 mm | 4U3005035090 | 90 mm  |
| 4U3005035048 | 48 mm | 4U3005035092 | 92 mm  |
| 4U3005035050 | 50 mm | 4U3005035094 | 94 mm  |
| 4U3005035052 | 52 mm | 4U3005035096 | 96 mm  |
| 4U3005035054 | 54 mm | 4U3005035098 | 98 mm  |
| 4U3005035055 | 55 mm | 4U3005035010 | 100 mm |
| 4U3005035056 | 56 mm | 4U3005035011 | 102 mm |
| 4U3005035058 | 58 mm | 4U3005035012 | 104 mm |
| 4U3005035060 | 60 mm | 4U3005035013 | 106 mm |
| 4U3005035062 | 62 mm | 4U3005035014 | 108 mm |
| 4U3005035064 | 64 mm | 4U3005035015 | 110 mm |

**LOCKING CANNULATED FULL  
THREAD SCREW Ø 6,5 MM**



**Code Size WASHER**

|              |        |
|--------------|--------|
| 4U3005006550 | 50 mm  |
| 4U3005006555 | 55 mm  |
| 4U3005006560 | 60 mm  |
| 4U3005006565 | 65 mm  |
| 4U3005006570 | 70 mm  |
| 4U3005006575 | 75 mm  |
| 4U3005006580 | 80 mm  |
| 4U3005006585 | 85 mm  |
| 4U3005006590 | 90 mm  |
| 4U3005006595 | 95 mm  |
| 4U3005006100 | 100 mm |
| 4U3005006105 | 105 mm |
| 4U3005006110 | 110 mm |
| 4U3005006115 | 115 mm |
| 4U3005006120 | 120 mm |



**Code Size**

|              |              |
|--------------|--------------|
| 4U4007000100 | S(3,5)       |
| 4U4007000101 | M(4,5)       |
| 4U4007000102 | L(6,5&7,3)   |
| 4U4007000103 | XL(6,5&7,3)  |
| 4U4007000104 | XXL(6,5&7,3) |

## HEADLESS CANNULATED FULL THREAD SCREW Ø 2,5 MM



| Code         | Size  |
|--------------|-------|
| 4U3009032008 | 8 mm  |
| 4U3009032012 | 10 mm |
| 4U3009032014 | 12 mm |
| 4U3009032016 | 14 mm |
| 4U3009032018 | 16 mm |
| 4U3009032020 | 18 mm |
| 4U3009032022 | 20 mm |
| 4U3009032024 | 22 mm |
| 4U3009032026 | 24 mm |
| 4U3009032028 | 26 mm |
| 4U3009032030 | 28 mm |
| 4U3009032010 | 30 mm |

## HEADLESS CANNULATED FULL THREAD SCREW Ø 3,5 MM



| Code         | Size  |
|--------------|-------|
| 4U3009032708 | 12 mm |
| 4U3009032710 | 14 mm |
| 4U3009032712 | 16 mm |
| 4U3009032714 | 18 mm |
| 4U3009032716 | 20 mm |
| 4U3009032718 | 22 mm |
| 4U3009032720 | 24 mm |
| 4U3009032722 | 26 mm |
| 4U3009032724 | 28 mm |
| 4U3009032726 | 30 mm |
| 4U3009032728 | 32 mm |
| 4U3009032730 | 34 mm |
| 4U3009033512 | 36 mm |
| 4U3009033514 | 38 mm |
| 4U3009033516 | 40 mm |
| 4U3009033518 | 42 mm |

## HEADLESS CANNULATED FULL THREAD SCREW Ø 4,0 MM




| Code         | Size  |
|--------------|-------|
| 4U3009033520 | 20 mm |
| 4U3009033522 | 22 mm |
| 4U3009033524 | 24 mm |
| 4U3009033526 | 26 mm |
| 4U3009033528 | 28 mm |
| 4U3009033530 | 30 mm |
| 4U3009033532 | 32 mm |
| 4U3009033534 | 34 mm |
| 4U3009033536 | 36 mm |
| 4U3009033538 | 38 mm |
| 4U3009033540 | 40 mm |
| 4U3009033542 | 45 mm |
| 4U3009033544 | 50 mm |
| 4U3009033546 | 55 mm |
| 4U3009033548 | 60 mm |
| 4U3009033550 | 65 mm |
| 4U3009034070 | 70 mm |

## HEADLESS CANNULATED FULL THREAD SCREW Ø 5,0 MM



| Code         | Size  | Code         | Size   |
|--------------|-------|--------------|--------|
| 4U3009035020 | 20 mm | 4U3009035065 | 65 mm  |
| 4U3009035025 | 25 mm | 4U3009035070 | 70 mm  |
| 4U3009035030 | 30 mm | 4U3009035075 | 75 mm  |
| 4U3009035035 | 35 mm | 4U3009035080 | 80 mm  |
| 4U3009035040 | 40 mm | 4U3009035085 | 85 mm  |
| 4U3009035045 | 45 mm | 4U3009035090 | 90 mm  |
| 4U3009035050 | 50 mm | 4U3009035095 | 95 mm  |
| 4U3009035055 | 55 mm | 4U3009035100 | 100 mm |
| 4U3009035060 | 60 mm |              |        |

## HERBERT SCREW Ø 5,0 MM



## HERBERT SCREW Ø 2,5 MM



## HERBERT SCREW Ø 3,0 MM



| Code         | Size  | Code         | Size  |
|--------------|-------|--------------|-------|
| 4U3009042508 | 8 mm  | 4U3009043008 | 8 mm  |
| 4U3009042509 | 9 mm  | 4U3009043010 | 10 mm |
| 4U3009042510 | 10 mm | 4U3009043012 | 12 mm |
| 4U3009042512 | 12 mm | 4U3009043014 | 14 mm |
| 4U3009042514 | 14 mm | 4U3009043016 | 16 mm |
| 4U3009042516 | 16 mm | 4U3009043018 | 18 mm |
| 4U3009042518 | 18 mm | 4U3009043020 | 20 mm |
| 4U3009042520 | 20 mm | 4U3009043022 | 22 mm |
| 4U3009042522 | 22 mm | 4U3009043024 | 24 mm |
| 4U3009042524 | 24 mm | 4U3009043026 | 26 mm |
| 4U3009042526 | 26 mm | 4U3009043028 | 28 mm |
| 4U3009042528 | 28 mm | 4U3009043030 | 30 mm |
| 4U3009042530 | 30 mm | 4U3009043032 | 32 mm |

## Code Size

|              |       |
|--------------|-------|
| 4U3009045020 | 20 mm |
| 4U3009045022 | 22 mm |
| 4U3009045024 | 24 mm |
| 4U3009045026 | 26 mm |
| 4U3009045028 | 28 mm |
| 4U3009045030 | 30 mm |
| 4U3009045032 | 32 mm |
| 4U3009045034 | 34 mm |
| 4U3009045036 | 36 mm |
| 4U3009045038 | 38 mm |
| 4U3009045040 | 40 mm |
| 4U3009045045 | 45 mm |
| 4U3009045050 | 50 mm |
| 4U3009045055 | 55 mm |
| 4U3009045060 | 60 mm |
| 4U3009045065 | 65 mm |
| 4U3009045070 | 70 mm |

### CABLE LOCKING SCREW



| Code         | Size |
|--------------|------|
| 4U3007035002 | 5 mm |

### CAP SCREW



| Code         | Size   |
|--------------|--------|
| 4U3007032500 | L(2,5) |
| 4U3007032501 | S(2,5) |
| 4U3007033500 | L(3,5) |
| 4U3007033501 | S(3,5) |
| 4U3007035000 | L(5,0) |
| 4U3007035001 | S(5,0) |

# 1,5 MM HAND / FOOT MINI LOCKING PLATE SERIES

## 1,5 Mini Locking Straight Plate

| Code         | Holes | Side | Thickness |
|--------------|-------|------|-----------|
| 4U1051001501 | 4     | -    | 1,2       |
| 4U1051001502 | 5     | -    | 1,2       |
| 4U1051001503 | 6     | -    | 1,2       |
| 4U1051001504 | 8     | -    | 1,2       |




## 1,5 Mini Locking T Plate Head 2 Holes

| Code         | Holes | Side | Thickness |
|--------------|-------|------|-----------|
| 4U1051001505 | 4     | -    | 1,2       |
| 4U1051001506 | 8     | -    | 1,2       |



## 1,5 Mini Locking H Plate

| Code         | Holes | Side | Thickness |
|--------------|-------|------|-----------|
| 4U1052001501 | -     | -    | 1,2       |



### Specification

Mini Plate Set  
Hand/Foot plate  
Intended Use : 1.5 mm mini locking plates are intended for temporary fixation, correction or stabilization of bones in anatomical region of the hand.

## 1,5 Mini Locking T Plate Head 3 Holes

| Code         | Holes | Side | Thickness |
|--------------|-------|------|-----------|
| 4U1051001507 | 4     | -    | 1,2       |
| 4U1051001508 | 8     | -    | 1,2       |



### Indications

Fracture fixation of the phalanges and metacarpals  
Osteotomies  
Arthrodesis  
Replantations and reconstructions of phalanges and metacarpals, particularly in osteopenic bone

## 1,5 Mini Locking T Plate Head 4 Holes

| Code         | Holes | Side | Thickness |
|--------------|-------|------|-----------|
| 4U1051002407 | 6     | -    | 1,2       |



## 1,5 Mini Locking Y Plate

|              | Holes | Side | Thickness |
|--------------|-------|------|-----------|
| 4U1051001509 | 4     | -    | 1,2       |
| 4U1051001510 | 8     | -    | 1,2       |



## 1,5 Mini Locking Condylar Plate

| Code         | Holes | Side | Thickness |
|--------------|-------|------|-----------|
| 4U1051001511 | 4     | -    | 1,2       |
| 4U1051001512 | 8     | -    | 1,2       |

### Specification

Mini Plate Set  
Hand/Foot plate  
Intended Use : 1.5 mm mini locking plates are intended for temporary fixation, correction or stabilization of bones in anatomical region of the hand.



### Indications

Fracture fixation of the phalanges and metacarpals  
Osteotomies  
Arthrodesis  
Replantations and reconstructions of phalanges and metacarpals, particularly in osteopenic bone

# 1,5 MM HAND / FOOT MINI LOCKING PLATE SERIES

## 1,5 Mini Non-Locking Strut Plate

| Code         | Holes | Side | Thickness |
|--------------|-------|------|-----------|
| 4U1051021513 | 8     | L    | 1,0       |
| 4U1051011514 | 8     | R    | 1,0       |

R



## 1,5 Mini Non-Locking Straight Plate

| Code         | Holes | Side | Thickness |
|--------------|-------|------|-----------|
| 4U1052001502 | 14    | -    | 1,0       |



## 1.5 Mini Non-locking Mallet Finger Plate

| Code         | Holes | Side | Thickness |
|--------------|-------|------|-----------|
| 4U1052001500 | 2     | -    | 1,0       |



Specification

Mini Plate Set  
Hand/Foot plate

Intended Use : 1.5 mm mini locking plates are intended for temporary fixation, correction or stabilization of bones in anatomical region of the hand.

## 1,5 Mini Non-Locking T Plate Head 4 Holes

| Code         | Holes | Side | Thickness |
|--------------|-------|------|-----------|
| 4U1051002406 | 6     | -    | 1,0       |



Indications

Fracture fixation of the phalanges and metacarpals  
Osteotomies  
Arthrodesis  
Replantations and reconstructions of phalanges and metacarpals, particularly in osteopenic bone

## LOCKING SELF TAPPING SCREW Ø1,5 MM



| Code         | Size  |
|--------------|-------|
| 4U3001001506 | 6 mm  |
| 4U3001001507 | 7 mm  |
| 4U3001001508 | 8 mm  |
| 4U3001001509 | 9 mm  |
| 4U3001001510 | 10 mm |
| 4U3001001511 | 11 mm |
| 4U3001001512 | 12 mm |
| 4U3001001513 | 13 mm |
| 4U3001001514 | 14 mm |
| 4U3001001515 | 15 mm |
| 4U3001001516 | 16 mm |
| 4U3001001517 | 17 mm |
| 4U3001001518 | 18 mm |
| 4U3001001519 | 19 mm |
| 4U3001001520 | 20 mm |
| 4U3001001522 | 22 mm |
| 4U3001001524 | 24 mm |
| 4U3001001526 | 26 mm |
| 4U3001001528 | 28 mm |
| 4U3001001530 | 30 mm |

## CORTICAL SELF TAPPING SCREW Ø1,5 MM



| Code         | Size  |
|--------------|-------|
| 4U3002001506 | 6 mm  |
| 4U3002001507 | 7 mm  |
| 4U3002001508 | 8 mm  |
| 4U3002001509 | 9 mm  |
| 4U3002001510 | 10 mm |
| 4U3002001511 | 11 mm |
| 4U3002001512 | 12 mm |
| 4U3002001513 | 13 mm |
| 4U3002001514 | 14 mm |
| 4U3002001515 | 15 mm |
| 4U3002001516 | 16 mm |
| 4U3002001517 | 17 mm |
| 4U3002001518 | 18 mm |
| 4U3002001519 | 19 mm |
| 4U3002001520 | 20 mm |
| 4U3002001522 | 22 mm |
| 4U3002001524 | 24 mm |
| 4U3002001526 | 26 mm |
| 4U3002001528 | 28 mm |
| 4U3002001530 | 30 mm |

# 2,0 MM HAND / FOOT MINI LOCKING PLATE SERIES

## 2.0 Mini Locking Straight Plate

| Code         | Holes | Side | Thickness |
|--------------|-------|------|-----------|
| 4U1051002000 | 4     | -    | 1,5       |
| 4U1051002001 | 5     | -    | 1,5       |
| 4U1051002002 | 6     | -    | 1,5       |
| 4U1051002003 | 7     | -    | 1,5       |
| 4U1051002004 | 8     | -    | 1,5       |



## 2,0 Mini Locking T Plate Head 2 Holes

| Code         | Holes | Side | Thickness |
|--------------|-------|------|-----------|
| 4U1051002007 | 3     | -    | 1,2       |
| 4U1051002008 | 5     | -    | 1,2       |
| 4U1051002009 | 7     | -    | 1,2       |



## 2,0 Mini Locking Condylar Plate

| Code         | Holes | Side | Thickness |
|--------------|-------|------|-----------|
| 4U1051002013 | 3     | -    | 1,2       |
| 4U1051002014 | 5     | -    | 1,2       |
| 4U1051002015 | 7     | -    | 1,2       |



### Specification

Mini Plate Set

Hand/Foot plate

Intended Use : 2.0 mm mini locking plates are intended for temporary fixation, correction or stabilization of bones in anatomical region of the hand.

## 2,0 Mini Locking Y Plate Head 3 Holes

| Code         | Holes | Side | Thickness |
|--------------|-------|------|-----------|
| 4U1051002010 | 3     | -    | 1,2       |
| 4U1051002011 | 5     | -    | 1,2       |
| 4U1051002012 | 7     | -    | 1,2       |



### Indications

Fracture fixation of the phalanges and metacarpals

Osteotomies

Arthrodesis

Replantations and reconstructions of phalanges and metacarpals, particularly in osteopenic bone

## 2.0 Mini Adaptation T Plate

| Code         | Holes | Side | Thickness |
|--------------|-------|------|-----------|
| 4U1051002016 | 3     | - -  | 1,2       |
| 4U1051002017 | 5     | -    | 1,2       |
| 4U1051002018 | 7     | -    | 1,2       |



## 2.0 Mini Adaptation Straight Plate

| Code         | Holes | Side | Thickness |
|--------------|-------|------|-----------|
| 4U1051002005 | 6     | -    | 1,2       |
| 4U1051002006 | 12    | -    | 1,2       |



## 2,0 MM HAND / FOOT MINI LOCKING PLATE SERIES

### LOCKING SELF TAPPING SCREW Ø2,0 MM



#### Code Size

|              |       |
|--------------|-------|
| 4U3001002006 | 6 mm  |
| 4U3001002007 | 7 mm  |
| 4U3001002008 | 8 mm  |
| 4U3001002009 | 9 mm  |
| 4U3001002010 | 10 mm |
| 4U3001002011 | 11 mm |
| 4U3001002012 | 12 mm |
| 4U3001002013 | 13 mm |
| 4U3001002014 | 14 mm |
| 4U3001002015 | 15 mm |
| 4U3001002016 | 16 mm |
| 4U3001002017 | 17 mm |
| 4U3001002018 | 18 mm |
| 4U3001002019 | 19 mm |
| 4U3001002020 | 20 mm |
| 4U3001002022 | 22 mm |
| 4U3001002024 | 24 mm |
| 4U3001002026 | 26 mm |
| 4U3001002028 | 28 mm |
| 4U3001002030 | 30 mm |

### CORTICAL SELF TAPPING SCREW Ø2,0 MM



#### Code Size

|              |       |
|--------------|-------|
| 4U3002002006 | 6 mm  |
| 4U3002002007 | 7 mm  |
| 4U3002002008 | 8 mm  |
| 4U3002002009 | 9 mm  |
| 4U3002002010 | 10 mm |
| 4U3002002011 | 11 mm |
| 4U3002002012 | 12 mm |
| 4U3002002013 | 13 mm |
| 4U3002002014 | 14 mm |
| 4U3002002015 | 15 mm |
| 4U3002002016 | 16 mm |
| 4U3002002017 | 17 mm |
| 4U3002002018 | 18 mm |
| 4U3002002019 | 19 mm |
| 4U3002002020 | 20 mm |
| 4U3002002022 | 22 mm |
| 4U3002002024 | 24 mm |
| 4U3002002026 | 26 mm |
| 4U3002002028 | 28 mm |
| 4U3002002030 | 30 mm |

## 2,5 MM HAND / FOOT MINI LOCKING PLATE SERIES

### 2,5 Mini Locking H Plate

#### Code Holes Side Thickness

|              |   |   |     |
|--------------|---|---|-----|
| 4U1051002517 | 5 | - | 1,5 |
| 4U1051002518 | 8 | - | 1,5 |



### 2,5 Mini Locking T Plate

#### Code Holes Side Thickness

|              |   |   |     |
|--------------|---|---|-----|
| 4U1051002519 | 3 | - | 1,3 |
| 4U1051002520 | 5 | - | 1,3 |
| 4U1051002521 | 7 | - | 1,3 |



# 2,5 MM HAND / FOOT MINI LOCKING PLATE SERIES

## 2.5 Mini Locking Straight Plate

| Code         | Holes | Side | Thickness |
|--------------|-------|------|-----------|
| 4U1051002500 | 4     | -    | 2,5       |
| 4U1051002501 | 5     | -    | 2,5       |
| 4U1051002502 | 6     | -    | 2,5       |
| 4U1051002503 | 7     | -    | 2,5       |
| 4U1051002504 | 8     | -    | 2,5       |
| 4U1051002505 | 10    | -    | 2,5       |



### Specification

Mini Plate Set  
Hand/Foot plate  
Intended Use : 2.5 mm mini locking plates are intended for temporary fixation, correction or stabilization of bones in anatomical region of the hand.

## 2,5 Mini Locking T Plate Head 2 Holes

| Code         | Holes | Side | Thickness |
|--------------|-------|------|-----------|
| 4U1051002408 | 3     | -    | 1,3       |
| 4U1051002409 | 5     | -    | 1,3       |
| 4U1051002410 | 7     | -    | 1,3       |



### Indications

Fracture fixation of the phalanges and metacarpals  
Osteotomies  
Arthrodesis  
Replantations and reconstructions of phalanges and metacarpals, particularly in osteopenic bone

## 2,5 Mini Locking Straight Plate

|              | Holes | Side | Thickness |
|--------------|-------|------|-----------|
| 4U1051002400 | 4     | -    | 2,0       |
| 4U1051002401 | 5     | -    | 2,0       |
| 4U1051002402 | 6     | -    | 2,0       |
| 4U1051002403 | 7     | -    | 2,0       |
| 4U1051002404 | 8     | -    | 2,0       |
| 4U1051002405 | 10    | -    | 2,0       |



## 2,5 Mini Locking T Plate Head 3 Holes

| Code         | Holes | Side | Thickness |
|--------------|-------|------|-----------|
| 4U1051002411 | 3     | -    | 1,5       |
| 4U1051002412 | 5     | -    | 1,5       |
| 4U1051002413 | 7     | -    | 1,5       |



## 2,5 Mini Locking L Plate

| Code         | Holes | Side | Thickness |
|--------------|-------|------|-----------|
| 4U1051002421 | 3     | R    | 1,5       |
| 4U1051002511 | 4     | R    | 1,5       |
| 4U1051002512 | 3     | L    | 1,5       |
| 4U1051002513 | 4     | L    | 1,5       |

R



### Specification

Mini Plate Set  
Hand/Foot plate  
Intended Use : 2.5 mm mini locking plates are intended for temporary fixation, correction or stabilization of bones in anatomical region of the hand.

## 2,5 Mini Adaptation Locking Straight Plate

| Code         | Holes | Side | Thickness |
|--------------|-------|------|-----------|
| 4U1051002506 | 6     | -    | 1,6       |
| 4U1051002507 | 12    | -    | 1,6       |



### Indications

Fracture fixation of the phalanges and metacarpals  
Osteotomies  
Arthrodesis  
Replantations and reconstructions of phalanges and metacarpals, particularly in osteopenic bone

## 2,5 MM HAND / FOOT MINI LOCKING PLATE SERIES

### 2,5 Mini Locking L Plate Oblique

| Code         | Holes | Side | Thickness |
|--------------|-------|------|-----------|
| 4U1051002417 | 3     | R    | 1,5       |
| 4U1051002418 | 4     | R    | 1,5       |
| 4U1051002419 | 3     | L    | 1,5       |
| 4U1051002420 | 4     | L    | 1,5       |

R



### 2,5 Mini Locking Condylar Plate

| Code         | Holes | Side | Thickness |
|--------------|-------|------|-----------|
| 4U1051002514 | 3     | -    | 1,3       |
| 4U1051002515 | 5     | -    | 1,3       |
| 4U1051002516 | 7     | -    | 1,3       |



### 2,5 Foot Locking X Plate

#### Code Holes Size Thickness

|              |   |   |     |
|--------------|---|---|-----|
| 4U1051002414 | 4 | S | 1,5 |
| 4U1051002415 | 4 | M | 1,5 |
| 4U1051002416 | 4 | L | 1,5 |



#### Specification

Mini Plate Set

Hand/Foot plate

Intended Use : 2.5 mm mini locking plates are intended for temporary fixation, correction or stabilization of bones in anatomical region of the hand.

#### Indications

Fracture fixation of the phalanges and metacarpals  
Osteotomies Arthrodesis Replantations and reconstructions of phalanges and metacarpals, particularly in osteopenic bone

### LOCKING SELF TAPPING SCREW Ø2,5 MM



### CORTICAL SELF TAPPING SCREW Ø2,5 MM



| Code         | Size  |
|--------------|-------|
| 4U3001002506 | 6 mm  |
| 4U3001002507 | 7 mm  |
| 4U3001002508 | 8 mm  |
| 4U3001002509 | 9 mm  |
| 4U3001002510 | 10 mm |
| 4U3001002511 | 11 mm |
| 4U3001002512 | 12 mm |
| 4U3001002513 | 13 mm |
| 4U3001002514 | 14 mm |
| 4U3001002515 | 15 mm |
| 4U3001002516 | 16 mm |
| 4U3001002517 | 17 mm |
| 4U3001002518 | 18 mm |
| 4U3001002519 | 19 mm |
| 4U3001002520 | 20 mm |
| 4U3001002522 | 22 mm |
| 4U3001002524 | 24 mm |
| 4U3001002526 | 26 mm |
| 4U3001002528 | 28 mm |
| 4U3001002530 | 30 mm |
| 4U3001002532 | 32 mm |
| 4U3001002534 | 34 mm |
| 4U3001002536 | 36 mm |
| 4U3001002538 | 38 mm |
| 4U3001002540 | 40 mm |

| Code         | Size  |
|--------------|-------|
| 4U3002002506 | 6 mm  |
| 4U3002002507 | 7 mm  |
| 4U3002002508 | 8 mm  |
| 4U3002002509 | 9 mm  |
| 4U3002002510 | 10 mm |
| 4U3002002511 | 11 mm |
| 4U3002002512 | 12 mm |
| 4U3002002513 | 13 mm |
| 4U3002002514 | 14 mm |
| 4U3002002515 | 15 mm |
| 4U3002002516 | 16 mm |
| 4U3002002517 | 17 mm |
| 4U3002002518 | 18 mm |
| 4U3002002519 | 19 mm |
| 4U3002002520 | 20 mm |
| 4U3002002522 | 22 mm |
| 4U3002002524 | 24 mm |
| 4U3002002526 | 26 mm |
| 4U3002002528 | 28 mm |
| 4U3002002530 | 30 mm |
| 4U3002002532 | 32 mm |
| 4U3002002534 | 34 mm |
| 4U3002002536 | 36 mm |
| 4U3002002538 | 38 mm |
| 4U3002002540 | 40 mm |

## 2,7 MM HAND / FOOT MINI LOCKING PLATE SERIES

### 2.7 Mini Locking Straight Plate

| Code         | Holes | Side | Thickness |
|--------------|-------|------|-----------|
| 4U1051002724 | 2     | -    | 2,5       |
| 4U1051002700 | 4     | -    | 2,5       |
| 4U1051002701 | 5     | -    | 2,5       |
| 4U1051002702 | 6     | -    | 2,5       |
| 4U1051002703 | 7     | -    | 2,5       |
| 4U1051002704 | 8     | -    | 2,5       |
| 4U1051002705 | 10    | -    | 2,5       |



### 2,7 Foot Locking X Plate

| Code         | Holes | Side | Thickness |
|--------------|-------|------|-----------|
| 4U1051002721 | 4     | S    | 1,8       |
| 4U1051002722 | 4     | M    | 1,8       |
| 4U1051002723 | 4     | L    | 1,8       |



### 2,7 Mini Locking H Plate

| Code         | Holes | Side | Thickness |
|--------------|-------|------|-----------|
| 4U1051002725 | 5     | -    | 1,8       |
| 4U1051002726 | 8     | -    | 1,8       |



### 2,7 Mini Locking Condylar Plate

| Code         | Holes | Side | Thickness |
|--------------|-------|------|-----------|
| 4U1051002718 | 3     | -    | 1,8       |
| 4U1051002719 | 5     | -    | 1,8       |
| 4U1051002720 | 7     | -    | 1,8       |



#### Specification

Mini Plate Set

Hand/Foot plate

Intended Use : 2.7 mm mini locking plates are intended for temporary fixation, correction or stabilization of bones in anatomical region of the hand.

#### Indications

Fracture fixation of the phalanges and metacarpals

Osteotomies

Arthrodesis

Replantations and reconstructions of phalanges and metacarpals, particularly in osteopenic bone

## 2,7 MM HAND / FOOT MINI LOCKING PLATE SERIES

### 2,7 Mini Locking L Plate

| Code         | Holes | Side | Thickness |
|--------------|-------|------|-----------|
| 4U1051012716 | 3     | R    | 1,8       |
| 4U1051012717 | 4     | R    | 1,8       |
| 4U1051022714 | 3     | L    | 1,8       |
| 4U1051022715 | 4     | L    | 1,8       |

R



### 2,7 Mini Locking T Plate

| Code         | Holes | Side | Thickness |
|--------------|-------|------|-----------|
| 4U1051002708 | 3     | -    | 1,8       |
| 4U1051002709 | 4     | -    | 1,8       |



### 2,7 Mini Locking L Plate Oblique

| Code         | Holes | Side | Thickness |
|--------------|-------|------|-----------|
| 4U1051012712 | 3     | R    | 1,8       |
| 4U1051012713 | 4     | R    | 1,8       |
| 4U1051022710 | 3     | L    | 1,8       |
| 4U1051022711 | 4     | L    | 1,8       |

R



#### Specification

Mini Plate Set

Hand/Foot plate

Intended Use : 2.7 mm mini locking plates are intended for temporary fixation, correction or stabilization of bones in anatomical region of the hand.

### 2,7 Mini Adaptation Locking Straight Plate

| Code         | Holes | Side | Thickness |
|--------------|-------|------|-----------|
| 4U1051002706 | 6     | -    | 1,8       |
| 4U1051002707 | 12    | -    | 1,8       |



#### Indications

Fracture fixation of the phalanges and metacarpals

Osteotomies

Arthrodesis

Replantations and reconstructions of phalanges and metacarpals, particularly in osteopenic bone

### LOCKING SELF TAPPING SCREW Ø2,7 MM



### CORTICAL SELF TAPPING SCREW Ø2,7 MM



| Code         | Size  | Code         | Size  | Code         | Size  |
|--------------|-------|--------------|-------|--------------|-------|
| 4U3001002710 | 10 mm | 4U3001002742 | 42 mm | 4U3002002710 | 10 mm |
| 4U3001002711 | 11 mm | 4U3001002744 | 44 mm | 4U3002002712 | 12 mm |
| 4U3001002712 | 12 mm | 4U3001002746 | 46 mm | 4U3002002714 | 14 mm |
| 4U3001002713 | 13 mm | 4U3001002748 | 48 mm | 4U3002002716 | 16 mm |
| 4U3001002714 | 14 mm | 4U3001002750 | 50 mm | 4U3002002718 | 18 mm |
| 4U3001002715 | 15 mm |              |       | 4U3002002720 | 20 mm |
| 4U3001002716 | 16 mm |              |       | 4U3002002722 | 22 mm |
| 4U3001002717 | 17 mm |              |       | 4U3002002724 | 24 mm |
| 4U3001002718 | 18 mm |              |       | 4U3002002726 | 26 mm |
| 4U3001002719 | 19 mm |              |       | 4U3002002728 | 28 mm |
| 4U3001002720 | 20 mm |              |       | 4U3002002730 | 30 mm |
| 4U3001002722 | 22 mm |              |       | 4U3002002732 | 32 mm |
| 4U3001002724 | 24 mm |              |       | 4U3002002734 | 34 mm |
| 4U3001002726 | 26 mm |              |       | 4U3002002736 | 36 mm |
| 4U3001002728 | 28 mm |              |       | 4U3002002738 | 38 mm |
| 4U3001002730 | 30 mm |              |       | 4U3002002740 | 40 mm |
| 4U3001002732 | 32 mm |              |       | 4U3002002742 | 42 mm |
| 4U3001002734 | 34 mm |              |       | 4U3002002744 | 44 mm |
| 4U3001002736 | 36 mm |              |       | 4U3002002746 | 46 mm |
| 4U3001002738 | 38 mm |              |       | 4U3002002748 | 48 mm |
| 4U3001002740 | 40 mm |              |       | 4U3002002750 | 50 mm |

METACARPAL / METATARSAL  
FUSION PLATE SERIES

| Code            | Product Description               | Angle                                      | Side | Thickness |                                                                                       |
|-----------------|-----------------------------------|--------------------------------------------|------|-----------|---------------------------------------------------------------------------------------|
| 4U1061202501    | Metacarpal Fusion Plate Narrow    | -                                          | -    | 2         |    |
| 4U1061202502    | Metacarpal Fusion Plate Broad     | -                                          | -    | 2         |    |
| 4U1091010200    | Mini Foot Metatarsal Fusion Plate | 0                                          | R    | 2         |    |
| 4U1091010201    | Mini Foot Metatarsal Fusion Plate | 5                                          | R    | 2         |    |
| 4U1091010202    | Mini Foot Metatarsal Fusion Plate | 15                                         | R    | 2         |                                                                                       |
| 4U1091020203    | Mini Foot Metatarsal Fusion Plate | 0                                          | L    | 2         |                                                                                       |
| 4U1091020204    | Mini Foot Metatarsal Fusion Plate | 5                                          | L    | 2         |  |
| 4U1091020205    | Mini Foot Metatarsal Fusion Plate | 15                                         | L    | 2         |                                                                                       |
| Specification   |                                   | Indications                                |      |           | R                                                                                     |
| Mini Plate Set  |                                   | Fractures of metatarsal                    |      |           |                                                                                       |
| Hand/Foot plate |                                   | Fractures of the tarsals                   |      |           |                                                                                       |
|                 |                                   | MTP 1 fusions                              |      |           |                                                                                       |
|                 |                                   | Osteotomies and arthrodesis of the tarsals |      |           |                                                                                       |

LOCKING SELF TAPPING SCREW  
Ø2,5 MM



CORTICAL SELF TAPPING SCREW  
Ø2,5 MM



| Code         | Size  | Code         | Size  | Code         | Size  | Code         | Size  |
|--------------|-------|--------------|-------|--------------|-------|--------------|-------|
| 4U3001002506 | 6 mm  | 4U3001002524 | 24 mm | 4U3002002506 | 6 mm  | 4U3002002524 | 24 mm |
| 4U3001002507 | 7 mm  | 4U3001002526 | 26 mm | 4U3002002507 | 7 mm  | 4U3002002526 | 26 mm |
| 4U3001002508 | 8 mm  | 4U3001002528 | 28 mm | 4U3002002508 | 8 mm  | 4U3002002528 | 28 mm |
| 4U3001002509 | 9 mm  | 4U3001002530 | 30 mm | 4U3002002509 | 9 mm  | 4U3002002530 | 30 mm |
| 4U3001002510 | 10 mm | 4U3001002532 | 32 mm | 4U3002002510 | 10 mm | 4U3002002532 | 32 mm |
| 4U3001002511 | 11 mm | 4U3001002534 | 34 mm | 4U3002002511 | 11 mm | 4U3002002534 | 34 mm |
| 4U3001002512 | 12 mm | 4U3001002536 | 36 mm | 4U3002002512 | 12 mm | 4U3002002536 | 36 mm |
| 4U3001002513 | 13 mm | 4U3001002538 | 38 mm | 4U3002002513 | 13 mm | 4U3002002538 | 38 mm |
| 4U3001002514 | 14 mm | 4U3001002540 | 40 mm | 4U3002002514 | 14 mm | 4U3002002540 | 40 mm |
| 4U3001002515 | 15 mm |              |       | 4U3002002515 | 15 mm |              |       |
| 4U3001002516 | 16 mm |              |       | 4U3002002516 | 16 mm |              |       |
| 4U3001002517 | 17 mm |              |       | 4U3002002517 | 17 mm |              |       |
| 4U3001002518 | 18 mm |              |       | 4U3002002518 | 18 mm |              |       |
| 4U3001002519 | 19 mm |              |       | 4U3002002519 | 19 mm |              |       |
| 4U3001002520 | 20 mm |              |       | 4U3002002520 | 20 mm |              |       |
| 4U3001002522 | 22 mm |              |       | 4U3002002522 | 22 mm |              |       |

# NAIL SYSTEM

## HUMERUS CANNULATED NAIL

|              |      |        | Humerus Locking Screw                                                               |                       |        |        |
|--------------|------|--------|-------------------------------------------------------------------------------------|-----------------------|--------|--------|
| Code         | Size | Length |  | Code                  | Size   | Length |
| 4U5024207160 | 7 mm | 160    |                                                                                     | 4U5026230020          | 5 mm   | 20     |
| 4U5024207185 | 7 mm | 180    |                                                                                     | 4U5026230025          | 5 mm   | 25     |
| 4U5024207200 | 7 mm | 200    |                                                                                     | 4U5026230030          | 5 mm   | 30     |
| 4U5024207215 | 7 mm | 220    |                                                                                     | 4U5026230035          | 5 mm   | 35     |
| 4U5024207245 | 7 mm | 240    |                                                                                     | 4U5026230040          | 5 mm   | 40     |
| 4U5024207260 | 7 mm | 260    |                                                                                     | 4U5026230045          | 5 mm   | 45     |
| 4U5024207275 | 7 mm | 280    |                                                                                     | 4U5026230050          | 5 mm   | 50     |
| 4U5024207305 | 7 mm | 300    |                                                                                     | 4U5026230055          | 5 mm   | 55     |
| 4U5024207320 | 7 mm | 320    |                                                                                     | 4U5026230060          | 5 mm   | 60     |
| 4U5024208160 | 8 mm | 160    |                                                                                     | 4U5026230065          | 5 mm   | 65     |
| 4U5024208185 | 8 mm | 180    |                                                                                     | 4U5026230070          | 5 mm   | 70     |
| 4U5024208200 | 8 mm | 200    |                                                                                     | 4U5026230075          | 5 mm   | 75     |
| 4U5024208215 | 8 mm | 220    |                                                                                     | 4U5026230080          | 5 mm   | 80     |
| 4U5024208245 | 8 mm | 240    |                                                                                     | Humerus Locking Screw |        |        |
| 4U5024208260 | 8 mm | 260    |                                                                                     | Code                  | Size   | Length |
| 4U5024208275 | 8 mm | 280    |                                                                                     | 4U5026235020          | 3,5 mm | 20     |
| 4U5024208305 | 8 mm | 300    |                                                                                     | 4U5026235025          | 3,5 mm | 25     |
| 4U5024208320 | 8 mm | 320    |                                                                                     | 4U5026235030          | 3,5 mm | 30     |
| 4U5024209160 | 9 mm | 160    |                                                                                     | 4U5026235035          | 3,5 mm | 35     |
| 4U5024209185 | 9 mm | 180    |                                                                                     | 4U5026235040          | 3,5 mm | 40     |
| 4U5024209200 | 9 mm | 200    |                                                                                     | 4U5026235045          | 3,5 mm | 45     |
| 4U5024209215 | 9 mm | 220    |                                                                                     | 4U5026235050          | 3,5 mm | 50     |
| 4U5024209245 | 9 mm | 240    |                                                                                     | 4U5026235055          | 3,5 mm | 55     |
| 4U5024209260 | 9 mm | 260    |                                                                                     | 4U5026235060          | 3,5 mm | 60     |
| 4U5024209275 | 9 mm | 280    |                                                                                     | 4U5026235065          | 3,5 mm | 65     |
| 4U5024209305 | 9 mm | 300    |                                                                                     | 4U5026235070          | 3,5 mm | 70     |
| 4U5024209320 | 9 mm | 320    |                                                                                     | 4U5026235075          | 3,5 mm | 75     |
|              |      |        |                                                                                     | 4U5026235080          | 3,5 mm | 80     |

### Indications

Stable or unstable fresh fractures, Pathologic fractures, Refractures, delayed unions and pseudarthrosis.

| Humerus End Cap |       |  |
|-----------------|-------|--|
| Code            | Size  |  |
| 4U5027030000    | 0     |  |
| 4U5026230005    | 5 mm  |  |
| 4U5026230010    | 10 mm |  |



# NAIL SYSTEM

## PFNA NAIL

| Code         | Size  | Side | Length |
|--------------|-------|------|--------|
| 4U5064109170 | 9 mm  | -    | 180    |
| 4U5064109200 | 9 mm  | -    | 200    |
| 4U5064109240 | 9 mm  | -    | 240    |
| 4U5064110170 | 10 mm | -    | 180    |
| 4U5064110200 | 10 mm | -    | 200    |
| 4U5064110240 | 10 mm | -    | 240    |
| 4U5064111170 | 11 mm | -    | 180    |
| 4U5064111200 | 11 mm | -    | 200    |
| 4U5064111240 | 11 mm | -    | 240    |
| 4U5064113170 | 13 mm | -    | 180    |
| 4U5064113200 | 13 mm | -    | 200    |
| 4U5064113240 | 13 mm | -    | 240    |
|              |       |      |        |
| 4U5064110300 | 10 mm | R    | 300    |
| 4U5064110320 | 10 mm | R    | 320    |
| 4U5064110340 | 10 mm | R    | 340    |
| 4U5064110360 | 10 mm | R    | 360    |
| 4U5064110380 | 10 mm | R    | 380    |
| 4U5064110400 | 10 mm | R    | 400    |
| 4U5064210300 | 10 mm | L    | 300    |
| 4U5064210320 | 10 mm | L    | 320    |
| 4U5064210340 | 10 mm | L    | 340    |
| 4U5064210360 | 10 mm | L    | 360    |
| 4U5064210380 | 10 mm | L    | 380    |
| 4U5064210400 | 10 mm | L    | 400    |
|              |       |      |        |
| 4U5064211300 | 11 mm | R    | 300    |
| 4U5064211320 | 11 mm | R    | 320    |
| 4U5064211340 | 11 mm | R    | 340    |
| 4U5064211360 | 11 mm | R    | 360    |
| 4U5064211380 | 11 mm | R    | 380    |
| 4U5064211400 | 11 mm | R    | 400    |
| 4U5064311300 | 11 mm | L    | 300    |
| 4U5064311320 | 11 mm | L    | 320    |
| 4U5064311340 | 11 mm | L    | 340    |
| 4U5064311360 | 11 mm | L    | 360    |
| 4U5064311380 | 11 mm | L    | 380    |
| 4U5064311400 | 11 mm | L    | 400    |



| Code         | Size  | Side | Length |
|--------------|-------|------|--------|
| 4U5064112300 | 12 mm | R    | 300    |
| 4U5064112320 | 12 mm | R    | 320    |
| 4U5064112340 | 12 mm | R    | 340    |
| 4U5064112360 | 12 mm | R    | 360    |
| 4U5064112380 | 12 mm | R    | 380    |
| 4U5064112400 | 12 mm | R    | 400    |
| 4U5064212300 | 12 mm | L    | 300    |
| 4U5064212320 | 12 mm | L    | 320    |
| 4U5064212340 | 12 mm | L    | 340    |
| 4U5064212360 | 12 mm | L    | 360    |
| 4U5064212380 | 12 mm | L    | 380    |
| 4U5064212400 | 12 mm | L    | 400    |
|              |       |      |        |
| 4U5064313300 | 13 mm | R    | 300    |
| 4U5064313320 | 13 mm | R    | 320    |
| 4U5064313340 | 13 mm | R    | 340    |
| 4U5064313360 | 13 mm | R    | 360    |
| 4U5064313380 | 13 mm | R    | 380    |
| 4U5064313400 | 13 mm | R    | 400    |
| 4U5064213300 | 13 mm | L    | 300    |
| 4U5064213320 | 13 mm | L    | 320    |
| 4U5064213340 | 13 mm | L    | 340    |
| 4U5064213360 | 13 mm | L    | 360    |
| 4U5064213380 | 13 mm | L    | 380    |
| 4U5064213400 | 13 mm | L    | 400    |

### Indications

Pertrochanteric fractures, Intertrochanteric fractures, High subtrochanteric fractures, Low subtrochanteric fractures, Femoral shaft fractures, Isolated or combined medial femoral neck fractures.

# NAIL SYSTEM

## PFNA NAIL

### PFNA Compression Locking Screw

| Code         | Length |
|--------------|--------|
| 4U5066030065 | 65     |
| 4U5066030070 | 70     |
| 4U5066030075 | 75     |
| 4U5066030080 | 80     |
| 4U5066030085 | 85     |
| 4U5066030090 | 90     |
| 4U5066030095 | 95     |
| 4U5066030100 | 100    |
| 4U5066030105 | 105    |
| 4U5066030110 | 110    |
| 4U5066030115 | 115    |
| 4U5066030120 | 120    |



### PFNA Anti-Rotation Locking Screw

| Code         | Length |
|--------------|--------|
| 4U7938140065 | 65     |
| 4U7938140070 | 70     |
| 4U7938140075 | 75     |
| 4U7938140080 | 80     |
| 4U7938140085 | 85     |
| 4U7938140090 | 90     |
| 4U7938140095 | 95     |
| 4U7938140100 | 100    |
| 4U7938140105 | 105    |
| 4U7938140110 | 110    |
| 4U7938140115 | 115    |
| 4U7938140120 | 120    |



### PFNA Nail End Cap Insert Locking Screw

| Code         | Size |
|--------------|------|
| 4U5007030025 | -    |



### PFNA Locking Screw

| Code         | Size | Length |
|--------------|------|--------|
| 4U5006230020 | 5 mm | 20     |
| 4U5006230025 | 5 mm | 25     |
| 4U5006230030 | 5 mm | 30     |
| 4U5006230035 | 5 mm | 35     |
| 4U5006230040 | 5 mm | 40     |
| 4U5006230045 | 5 mm | 45     |
| 4U5006230050 | 5 mm | 50     |
| 4U5006230055 | 5 mm | 55     |
| 4U5006230060 | 5 mm | 60     |
| 4U5006230065 | 5 mm | 65     |
| 4U5006230070 | 5 mm | 70     |
| 4U5006230075 | 5 mm | 75     |
| 4U5006230080 | 5 mm | 80     |
| 4U5006230085 | 5 mm | 85     |
| 4U5006230090 | 5 mm | 90     |
| 4U5006230095 | 5 mm | 95     |
| 4U5006230100 | 5 mm | 100    |
| 4U5006230105 | 5 mm | 105    |
| 4U5006230110 | 5 mm | 110    |
| 4U5006230115 | 5 mm | 115    |
| 4U5006230120 | 5 mm | 120    |



### PFNA Nail End Cap

| Code         | Size  |
|--------------|-------|
| 4U5007030000 | 0     |
| 4U5007030005 | 5 mm  |
| 4U5007030010 | 10 mm |
| 4U5007030015 | 15 mm |
| 4U5007030020 | 20 mm |



## NAIL SYSTEM

### INTRAMEDULLARY CANNULATED LATERAL FEMUR NAIL

| Code         | Size   | Side | Length |                                                                                    | Code         | Size    | Side | Length |
|--------------|--------|------|--------|------------------------------------------------------------------------------------|--------------|---------|------|--------|
| 4U5064318300 | 9,5 mm | R    | 300    |  | 4U5064111300 | 11,5 mm | R    | 300    |
| 4U5064318320 | 9,5 mm | R    | 320    |                                                                                    | 4U5064111320 | 11,5 mm | R    | 320    |
| 4U5064318340 | 9,5 mm | R    | 340    |                                                                                    | 4U5064111340 | 11,5 mm | R    | 340    |
| 4U5064318360 | 9,5 mm | R    | 360    |                                                                                    | 4U5064111360 | 11,5 mm | R    | 360    |
| 4U5064318380 | 9,5 mm | R    | 380    |                                                                                    | 4U5064111380 | 11,5 mm | R    | 380    |
| 4U5064318400 | 9,5 mm | R    | 400    |                                                                                    | 4U5064111400 | 11,5 mm | R    | 400    |
| 4U5064318420 | 9,5 mm | R    | 420    |                                                                                    | 4U5064111420 | 11,5 mm | R    | 420    |
| 4U5064318440 | 9,5 mm | R    | 440    |                                                                                    | 4U5064111440 | 11,5 mm | R    | 440    |
| 4U5064328400 | 9,5 mm | L    | 300    |                                                                                    | 4U5064221300 | 11,5 mm | L    | 300    |
| 4U5064328420 | 9,5 mm | L    | 320    |                                                                                    | 4U5064221320 | 11,5 mm | L    | 320    |
| 4U5064328440 | 9,5 mm | L    | 340    |                                                                                    | 4U5064221340 | 11,5 mm | L    | 340    |
| 4U5064328460 | 9,5 mm | L    | 360    |                                                                                    | 4U5064221360 | 11,5 mm | L    | 360    |
| 4U5064328480 | 9,5 mm | L    | 380    |                                                                                    | 4U5064221380 | 11,5 mm | L    | 380    |
| 4U5064328500 | 9,5 mm | L    | 400    |                                                                                    | 4U5064221400 | 11,5 mm | L    | 400    |
| 4U5064328520 | 9,5 mm | L    | 420    |                                                                                    | 4U5064221420 | 11,5 mm | L    | 420    |
| 4U5064328540 | 9,5 mm | L    | 440    |                                                                                    | 4U5064221440 | 11,5 mm | L    | 440    |
| 4U5064118300 | 10 mm  | R    | 300    |                                                                                    | 4U5064113300 | 12 mm   | R    | 300    |
| 4U5064118320 | 10 mm  | R    | 320    |                                                                                    | 4U5064113320 | 12 mm   | R    | 320    |
| 4U5064118340 | 10 mm  | R    | 340    |                                                                                    | 4U5064113340 | 12 mm   | R    | 340    |
| 4U5064118360 | 10 mm  | R    | 360    |                                                                                    | 4U5064113360 | 12 mm   | R    | 360    |
| 4U5064118380 | 10 mm  | R    | 380    |                                                                                    | 4U5064113380 | 12 mm   | R    | 380    |
| 4U5064118400 | 10 mm  | R    | 400    |                                                                                    | 4U5064113400 | 12 mm   | R    | 400    |
| 4U5064110420 | 10 mm  | R    | 420    |                                                                                    | 4U5064113420 | 12 mm   | R    | 420    |
| 4U5064110440 | 10 mm  | R    | 440    |                                                                                    | 4U5064113440 | 12 mm   | R    | 440    |
| 4U5064220300 | 10 mm  | L    | 300    |                                                                                    | 4U5064223300 | 12 mm   | L    | 300    |
| 4U5064220320 | 10 mm  | L    | 320    |                                                                                    | 4U5064223320 | 12 mm   | L    | 320    |
| 4U5064220340 | 10 mm  | L    | 340    |                                                                                    | 4U5064223340 | 12 mm   | L    | 340    |
| 4U5064220360 | 10 mm  | L    | 360    |                                                                                    | 4U5064223360 | 12 mm   | L    | 360    |
| 4U5064220380 | 10 mm  | L    | 380    |                                                                                    | 4U5064223380 | 12 mm   | L    | 380    |
| 4U5064220400 | 10 mm  | L    | 400    |                                                                                    | 4U5064223400 | 12 mm   | L    | 400    |
| 4U5064220420 | 10 mm  | L    | 420    |                                                                                    | 4U5064223420 | 12 mm   | L    | 420    |
| 4U5064220440 | 10 mm  | L    | 440    |                                                                                    | 4U5064223440 | 12 mm   | L    | 440    |

#### Indications

Lateral Femoral Nail with standard locking is indicated for fractures in the femoral shaft except subtrochanteric fractures femoral shaft in case of combination with femoral neck fractures: combined with double ipsilateral fractures indicated for fractures in the subtrochanteric section.

## NAIL SYSTEM

### INTRAMEDULLARY CANNULATED LATERAL FEMUR NAIL

#### Lateral Femur Nail Cannulated Lag Screw

| Code         | Side Length |     |
|--------------|-------------|-----|
| 4U5066230065 | 6,5         | 65  |
| 4U5066230070 | 6,5         | 70  |
| 4U5066230075 | 6,5         | 75  |
| 4U5066230080 | 6,5         | 80  |
| 4U5066230085 | 6,5         | 85  |
| 4U5066230090 | 6,5         | 90  |
| 4U5066230095 | 6,5         | 95  |
| 4U5066230100 | 6,5         | 100 |
| 4U5066230105 | 6,5         | 105 |
| 4U5066230110 | 6,5         | 110 |
| 4U5066230115 | 6,5         | 115 |
| 4U5066230120 | 6,5         | 120 |



#### Lateral Femur Nail End Cap

| Code         | Side  |
|--------------|-------|
| 4U5007031000 | 0     |
| 4U5007031005 | 5 mm  |
| 4U5007031010 | 10 mm |
| 4U5007031015 | 15 mm |
| 4U5007031020 | 20 mm |



#### Lateral Femur Nail Locking Screw

| Code         | Size | Length |
|--------------|------|--------|
| 4U5064205020 | 5 mm | 20     |
| 4U5064205025 | 5 mm | 25     |
| 4U5064205028 | 5 mm | 28     |
| 4U5064205030 | 5 mm | 30     |
| 4U5064205032 | 5 mm | 32     |
| 4U5064205034 | 5 mm | 34     |
| 4U5064205035 | 5 mm | 35     |
| 4U5064205036 | 5 mm | 36     |
| 4U5064205038 | 5 mm | 38     |
| 4U5064205040 | 5 mm | 40     |
| 4U5064205042 | 5 mm | 42     |
| 4U5064205044 | 5 mm | 44     |
| 4U5064205045 | 5 mm | 45     |
| 4U5064205048 | 5 mm | 48     |
| 4U5064205050 | 5 mm | 50     |
| 4U5064205052 | 5 mm | 52     |
| 4U5064205054 | 5 mm | 54     |
| 4U5064205055 | 5 mm | 55     |
| 4U5064205056 | 5 mm | 56     |
| 4U5064205058 | 5 mm | 58     |
| 4U5064205060 | 5 mm | 60     |
| 4U5064205064 | 5 mm | 64     |
| 4U5064205065 | 5 mm | 65     |
| 4U5064205068 | 5 mm | 68     |
| 4U5064205070 | 5 mm | 70     |
| 4U5064205072 | 5 mm | 72     |
| 4U5064205074 | 5 mm | 74     |
| 4U5064205075 | 5 mm | 75     |
| 4U5064205076 | 5 mm | 76     |
| 4U5064205078 | 5 mm | 78     |
| 4U5064205080 | 5 mm | 80     |
| 4U5064205085 | 5 mm | 85     |
| 4U5064205090 | 5 mm | 90     |
| 4U5064205095 | 5 mm | 95     |
| 4U5064205100 | 5 mm | 100    |
| 4U5064205105 | 5 mm | 105    |
| 4U5064205110 | 5 mm | 110    |
| 4U5064205115 | 5 mm | 115    |
| 4U5064205120 | 5 mm | 120    |



# NAIL SYSTEM

## TIBIAL NAIL CANNULATED

| Code         | Size   | Length | Code         | Size  | Length |
|--------------|--------|--------|--------------|-------|--------|
| 4U5074208255 | 8,4 mm | 255    | 4U5074211255 | 11 mm | 255    |
| 4U5074208270 | 8,4 mm | 270    | 4U5074211270 | 11 mm | 270    |
| 4U5074208280 | 8,4 mm | 285    | 4U5074211280 | 11 mm | 285    |
| 4U5074208300 | 8,4 mm | 300    | 4U5074211300 | 11 mm | 300    |
| 4U5074208315 | 8,4 mm | 315    | 4U5074211320 | 11 mm | 315    |
| 4U5074208330 | 8,4 mm | 330    | 4U5074211330 | 11 mm | 330    |
| 4U5074208345 | 8,4 mm | 345    | 4U5074211340 | 11 mm | 345    |
| 4U5074208360 | 8,4 mm | 360    | 4U5074211360 | 11 mm | 360    |
| 4U5074208375 | 8,4 mm | 375    | 4U5074211380 | 11 mm | 375    |
| 4U5074208400 | 8,4 mm | 400    | 4U5074211400 | 11 mm | 400    |
| 4U5074209255 | 9 mm   | 255    | 4U5074212270 | 12 mm | 270    |
| 4U5074209270 | 9 mm   | 270    | 4U5074212280 | 12 mm | 285    |
| 4U5074209280 | 9 mm   | 285    | 4U5074212300 | 12 mm | 300    |
| 4U5074209300 | 9 mm   | 300    | 4U5074212320 | 12 mm | 315    |
| 4U5074209320 | 9 mm   | 315    | 4U5074212330 | 12 mm | 330    |
| 4U5074209330 | 9 mm   | 330    | 4U5074212340 | 12 mm | 345    |
| 4U5074209340 | 9 mm   | 345    | 4U5074212360 | 12 mm | 360    |
| 4U5074209360 | 9 mm   | 360    | 4U5074212380 | 12 mm | 375    |
| 4U5074209380 | 9 mm   | 375    | 4U5074212400 | 12 mm | 400    |
| 4U5074209400 | 9 mm   | 400    |              |       |        |
| 4U5074210255 | 10 mm  | 255    |              |       |        |
| 4U5074210270 | 10 mm  | 270    |              |       |        |
| 4U5074210280 | 10 mm  | 285    |              |       |        |
| 4U5074210300 | 10 mm  | 300    |              |       |        |
| 4U5074210320 | 10 mm  | 315    |              |       |        |
| 4U5074210330 | 10 mm  | 330    |              |       |        |
| 4U5074210340 | 10 mm  | 345    |              |       |        |
| 4U5074210360 | 10 mm  | 360    |              |       |        |
| 4U5074210380 | 10 mm  | 375    |              |       |        |
| 4U5074210400 | 10 mm  | 400    |              |       |        |



### Indications


Fractures in the tibial shaft as well as for metaphyseal and certain intraarticular fractures of the tibial head and the pilon tibial.

## NAIL SYSTEM

### TIBIAL NAIL CANNULATED

#### Tibial Nail Locking Screw

| Code         | Size   | Length |
|--------------|--------|--------|
| 4U5076235025 | 4,8 mm | 25     |
| 4U5076235030 | 4,8 mm | 30     |
| 4U5076235035 | 4,8 mm | 35     |
| 4U5076235040 | 4,8 mm | 40     |
| 4U5076235045 | 4,8 mm | 45     |
| 4U5076235050 | 4,8 mm | 50     |
| 4U5076235055 | 4,8 mm | 55     |
| 4U5076235060 | 4,8 mm | 60     |
| 4U5076235065 | 4,8 mm | 65     |
| 4U5076235070 | 4,8 mm | 70     |
| 4U5076235075 | 4,8 mm | 75     |
| 4U5076235080 | 4,8 mm | 80     |
| 4U5076235085 | 4,8 mm | 85     |
| 4U5076235090 | 4,8 mm | 90     |
| 4U5076235095 | 4,8 mm | 95     |
| 4U5076235010 | 4,8 mm | 100    |



#### Tibial Nail Locking Screw

| Code         | Size   | Length |
|--------------|--------|--------|
| 4U5076234520 | 4,5 mm | 20     |
| 4U5076234524 | 4,5 mm | 24     |
| 4U5076234525 | 4,5 mm | 25     |
| 4U5076234526 | 4,5 mm | 26     |
| 4U5076234528 | 4,5 mm | 28     |
| 4U5076234530 | 4,5 mm | 30     |
| 4U5076234532 | 4,5 mm | 32     |
| 4U5076234534 | 4,5 mm | 34     |
| 4U5076234535 | 4,5 mm | 35     |
| 4U5076234536 | 4,5 mm | 36     |
| 4U5076234538 | 4,5 mm | 38     |
| 4U5076234540 | 4,5 mm | 40     |
| 4U5076234542 | 4,5 mm | 42     |
| 4U5076234544 | 4,5 mm | 44     |
| 4U5076234545 | 4,5 mm | 45     |
| 4U5076234546 | 4,5 mm | 46     |
| 4U5076234548 | 4,5 mm | 48     |
| 4U5076234550 | 4,5 mm | 50     |
| 4U5076234554 | 4,5 mm | 54     |
| 4U5076234555 | 4,5 mm | 55     |
| 4U5076234560 | 4,5 mm | 60     |
| 4U5076234565 | 4,5 mm | 65     |
| 4U5076234570 | 4,5 mm | 70     |
| 4U5076234575 | 4,5 mm | 75     |
| 4U5076234580 | 4,5 mm | 80     |
| 4U5076234585 | 4,5 mm | 85     |
| 4U5076234590 | 4,5 mm | 90     |
| 4U5076234595 | 4,5 mm | 95     |
| 4U5076234510 | 4,5 mm | 100    |



#### Tibial Nail End

| Code         | Cap Size |
|--------------|----------|
| 4U5007032000 | 0        |
| 4U5007032005 | 5 mm     |
| 4U5007032010 | 10 mm    |
| 4U5007032015 | 15 mm    |
| 4U5007032020 | 20 mm    |



# NAIL SYSTEM

## ARTHRODESIS NAIL

| Code         | Size  | Length |
|--------------|-------|--------|
| 4U8501000150 | 10 mm | 150    |
| 4U8501000175 | 10 mm | 175    |
| 4U8501000200 | 10 mm | 200    |
| 4U8501000225 | 10 mm | 225    |
| 4U8501100150 | 11 mm | 150    |
| 4U8501100175 | 11 mm | 175    |
| 4U8501100200 | 11 mm | 200    |
| 4U8501100225 | 11 mm | 225    |

### Arthrodesis Locking Screw

| Code         | Size   | Length |
|--------------|--------|--------|
| 4U6026230020 | 5.0 mm | 20     |
| 4U6026230025 | 5.0 mm | 25     |
| 4U6026230030 | 5.0 mm | 30     |
| 4U6026230035 | 5.0 mm | 35     |
| 4U6026230040 | 5.0 mm | 40     |
| 4U6026230045 | 5.0 mm | 45     |
| 4U6026230050 | 5.0 mm | 50     |
| 4U6026230055 | 5.0 mm | 55     |
| 4U6026230060 | 5.0 mm | 60     |
| 4U6026230065 | 5.0 mm | 65     |
| 4U6026230070 | 5.0 mm | 70     |
| 4U6026230075 | 5.0 mm | 75     |
| 4U6026230080 | 5.0 mm | 80     |

### Arthrodesis Nail End Cap

| Code         | Size |
|--------------|------|
| 4U6026230085 | 0    |

#### Indications

Pertrochanteric fractures, Intertrochanteric fractures, High subtrochanteric fractures, Low subtrochanteric fractures, Femoral shaft fractures, Isolated or combined medial femoral neck fractures.

## FLEXIBLE INTRAMEDULLARY NAIL

| Code         | Size   | Length |
|--------------|--------|--------|
| 4U5032202040 | 2,0 mm | 400    |
| 4U5032202540 | 2,5 mm | 400    |
| 4U5032203040 | 3,0 mm | 400    |
| 4U5032203540 | 3,5 mm | 400    |
| 4U5032204040 | 4,0 mm | 400    |
| 4U5032204540 | 4,5 mm | 400    |
| 4U5032205040 | 5,0 mm | 400    |

### Flexible Nail End Cap

| Code         | Size |
|--------------|------|
| 4U5032202041 | 0    |
| 4U5037030000 | 0    |

#### Indications

The management of diaphyseal and certain metaphyseal/-epiphyseal fractures of long bones in children and young adults.



MEDICAL  
4U®

Turkey



[www.4umedical.com.tr](http://www.4umedical.com.tr)

# Keeps You Moving



**4U Mak. Medical  
Oto.San.Tic.Ltd.Şti.**

-  Turgut Ozal Mah. 2118.Cd.No:6/A  
Yenimahalle/Ankara-TURKIYE
-  +90 312 566 06 55
-  +90 532 677 42 32
-  [info@4umedical.com.tr](mailto:info@4umedical.com.tr)



[www.4umedical.com.tr](http://www.4umedical.com.tr)

